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ENGLISH

INDEX

1. Introduction	69
1.1 Areas of use/indications for Coyote T+	70
1.2 Contra indications	71
1.3 Quality and durability	71
1.4 The environment and waste disposal	71
1.5 Information for re-use	72
1.6 About this manual	72
1.7 Vital measures	73
2. Quick reference	76
3. Description	77
3.1 Information about Coyote T+	77
4. Features	78
4.1 Accessories	78
5. Assembling and adjusting	80
5.1 Unpacking	80
5.2 Initial adjustments to the user - overview	80
5.3 Wheels assembling	82
5.3.1 16" Drive wheels	82
5.3.2 Drive wheels with quick release	82
5.3.3 Correct frame angle	83
5.3.4 Main wheel quick release	83
5.3.5 Repositioning of the main wheel brackets for- and backwards	83
5.3.6 Height adjustments by repositioning wheels in the main wheels brackets	84
5.3.7 Anti-tips	86
5.3.8 Stepper pedal adjustment	87
5.3.9 Front fork	87
5.3.10 Castor position in front fork	87
5.3.11 Front castors	88
5.3.12 The brakes	88
5.3.13 User brake adjustment	89
5.3.14 Drum brake adjustment	89
5.4 Seat assembling - wheelchair seat height	90
5.4.1 Seat height	90
5.4.2 Seat plate - seat depth adjustments	90
5.5 Chair back assembling and positioning	93
5.6. Assessment and choice of back support gas spring	93
5.6.1 Change of recline gas spring	94
5.6.2 Change of tilt gas spring	94
5.7 Head support assembling and positioning	95
5.7.1 Head supports	95

5.8 Arm supports adjusting	96
5.8.1 Arm supports height	96
5.8.2 Arm supports depth	96
5.8.3 Arm supports handling	97
5.9 Leg supports	97
5.9.1 Leg support height and length adjustments	97
5.9.2 Universal leg supports	98
5.9.3 Angle adjustable leg supports	99
5.9.4 Leg support mounting and removal	101
5.9.5 Foot plate lock	101
5.10 Abduction block	101
6. Everyday use	102
6.1 Key words regarding tilt and recline of comfort wheelchairs	102
6.2 Decrease the possibility of sliding, shear and pressure sores	103
6.3 Operating recline: reclining the chair back	103
6.4 Operating tilt: tilting the seating unit	104
6.5 Practice the use, tilt and recline system	105
6.6 Transfer to and from the wheelchair	105
6.7 Using a patient hoist	105
6.8 Daily attendant handling	106
6.8.1 Push bow	106
6.8.2 Anti-tips	106
6.8.3 Tramping pedal	106
6.8.4 Brakes	107
6.8.5 Arm supports	107
6.8.6 Head support-mounting	108
6.8.7 Leg supports - mounting	109
7. Transport	110
7.1 Transport in car	110
7.2 Folding for transport	112
7.3 Transport in air-plane	112
7.4 Travelling on public transport	113
8. Manoeuvring	113
8.1 General techniques	113
8.2 Driving techniques - step up	115
8.3 Driving techniques - step down	115
8.4 Driving techniques - slope	116
8.5 Driving techniques - up and down stairs	116
8.6 Transfer	116
8.7 Point of balance	117
8.8 Lifting the wheelchair	118
9. Maintenance	118
9.1 Maintenance instructions	118
9.2 Cleaning and washing	119

9.3 Long term storing	119
10. Troubleshooting	120
11. Tests & Warranty	121
11.1 Tests	121
11.2 Warranty	121
11.3 Claim	121
11.4 Customized / individual adaptations	122
11.5 Combinations with other products	122
11.6 Service and repair	123
12. Measurements & weight	124

Dealer:
Serial number:
Date:
Stamp:

1. INTRODUCTION

Dear Customer,

You have made a great choice by purchasing a **Coyote T+** tilt-in-space wheelchair!

On the following pages, you will learn how to make the most of everything this product has to offer.

Thank you very much for your trust!

This wheelchair has been tested and certified in accordance with the requirements of FDA, CE, BSEN 12183 and Standards ABNT NBR ISO 7176-1 :2018, 7176-3 :2015, 7176-7 :2009, 7176-8 :2009, 7176-15 :2009.

Manufacturer:

Loh Enterprises, LLC dba Loh Medical
318 Davis St.
Clarks Summit, PA 18411
USA

To contact customer service, visit:

<https://www.lohmedical.com/en/contact>

Coyote T+ is a tilt and recline wheelchair for adolescent, adults and grown up users. The tests were carried out by TÜV SÜD Product Service GmbH in Germany.

In OneLoh we believe that wheelchairs should be chosen based on a thorough assessment focusing on the needs of the user and demands from environment. It is vital to be aware of the possibilities and limits the wheelchair offers.

Coyote T+ is a wheelchair designed for users with great need of comfort and adjustments. It has an adjustable seat and back angle, thus facilitating for the user change of position, mobilisation or posture correction (stabilization).

Coyote T+ has been crash tested at TASS International test centre in Helmond Netherlands, according to ISO 7176-19:2022

Max user weight is 136 kg., also when used as a seat in a car with transit system included. When mounting accessories such as power kit etc, the weight of the accessories must be subtracted from the max user weight. Specifications varies between countries.

Coyote T+ wheelchair is identified in the Global Medical Device Nomenclature GMDN 41620 -

Wheelchair attendant/occupant driven, rear wheels driven, non-collapsible.

1.1 Areas of use/indications

Coyote T+ is a multifunctional, tilt and recline comfort wheelchair for wheeled transport for partially or fully immobile adolescent and grown up persons with physical and/or mental disabilities. **Coyote T+** has a modular system that can be customized and adjusted according to the user's need.

It provides a means by which a disabled person who requires the possibility to vary the sitting position from activity to rest will profit from using tilt and recline functions. It is designed for users with the need for comfort and relief.

The combination between the seating system and the ergonomic solution in the frame construction offers many possibilities for adaptation and adjustments.

IMPORTANT BENEFITS of Coyote T+

The disabilities may have multiple causes. **Coyote T+** has an adjustable seat and back angle, thus facilitating for the user change of position, mobilization or posture correction (stabilization), wherever the following functional impediments with their multiple possible causes are present:

- Limited or lacking mobility
- Limited or lacking muscle power

- Limited movement range
- Lacking or limited trunk and body stability
- Hemiplegia
- Rheumatic-type disorders
- Craniocerebral injuries
- Cerebral Palsy and other neurological or geriatric disorders.

Before using **Coyote T+**, a seating assessment should be carried out by a trained professional.

If our standard solutions does not cover you needs, please contact our customer service for customized solutions.

1.2 Contra indications

Should not be used if:

- The user has strong involuntary movements. It is recommended the use of Dynamic System for these users.

1.3 Quality and durability

Coyote T+ wheelchair is tested at TÜV SÜD Product Service GmbH in Germany, following the European Standard DIN EN 12183:2014.

As manufacturer, OneLoh A.S evaluates the test to be equal to 5-6 years of normal use of the chair.

The disability of the user, the toughness of use as well as the level and regularity of maintainance done, foremost decides the durability of the wheelchair.

Thus, the durability will vary depending on these factors. With adequate maintenance, the lifetime of the wheelchair can be expected to exceed the 5 year warranty period by many years.

1.4 The environment and waste disposal

OneLoh and its suppliers wish to protect the environment.

This means:

- That we avoid using environmentally harmful substances and processes to the greatest extent possible.
- That OneLoh's products are ensured a long service life and a high degree of flexibility - to benefit the environment and economy.
- That all packaging can be recycled.
- That the wheelchair was designed to be separated into its components materials - to make recycling easier.

Contact your local recycling agent to get correct information how to handle in your area.

TEMPERATURE RANGE

Coyote T+ wheelchair is designed for temperature range of -10°C to +40°C

1.5 Information for re-use

All products from OneLoh are designed to require little maintenance and service over years. All products can be adapted for re-use by an authorised dealer. In order to guarantee performance and safety, OneLoh recommends the following tests prior to any re-use.

Please examine the following components for function, integrity etc. and replace parts if necessary:

- Wheelchair frame
- Front castors
- Solid rear wheels
- Hubs
- Brake function
- Directional stability of wheels
- Bearings: test for wear and lubrication.
- Leg supports
- Arm supports
- Recline / tilt function
- Push bar / handles
- Anti tip

For hygienic reasons: please replace the head support cover for a new user. Please also note the contents of Section 9.2 Cleaning and care.

1.6 About this manual

In order to avoid damages while using the **Coyote T+**, please read this manual carefully before starting to use the chair.

Symbol of forbidden actions.

No warranty can be claimed whenever these actions are implemented.

Symbol of warning.

Whenever this symbol is used, caution has to be taken.

Symbol for important information.

Symbol for useful tips.

Symbol for tools.

Symbol for:

Max safe slope for hand brake.

Symbol for: Max user weight. 136 kg.

Medical Device

Manufacturer - Name and address

Date of manufacturing

Serial number

Read user instruction

Please note that this manual is updated according to the year and date stated on each page.

User Manual on web -lohmedical.com

For enhanced readability (advantageous for users with visibility challenges) please find our user manual on our web page: www.lohmedical.com -manuals - user manual **Coyote T+**.

Latest user manual updates, product safety notes, addresses and other product information like recalls etc. will be published on our web page.

1.7 Vital measures

Coyote T+ is a multifunctional class B wheelchair designed for indoor and outdoor use.

Min. dimensions in table refer to seat width 350 mm. Max dim refers to seat width 500 mm.

Specifications varies between countries.

Total weight: 34 kg.
(450 mm width chair)


Seat width

350, 380, 400, 430, 450, 500 mm

(8mm less between arm support pads):

Overall width:

690mm for chair width 500mm by main wheels 16",

Seat depth

From back to front of seat plate:

350, 400, 425, 450, 475, 500 mm


Seat height:

435 mm*

From floor to top seat plate dependent on wheel sizes and hole position.

*With same wheel sizes the seating height can be adjusted in the wheel brackets.

Backrest height:

515 mm*

* Using the back rest extender gives up to 610 mm backrest height.

For Coyote T+ with seat width 500 mm, the overall width exceeds recommended 700 mm.

Specification Coyote T+ - seat width 450 mm	min.	max.
Overall length with leg rest – horizontal seat. (Seat + back tilted max forward. Leg supports vertical)	1010 mm (946 mm)	1150 mm
Overall width (230 mm + seat width)**	680 mm	-
Folded length	780 mm	780 mm
Folded width (removed wheels)	640 mm	-
Folded height (removed wheels, cushions, arm, head + leg support)	578 mm	648 mm
Total mass	34 kg	34 kg
Mass heaviest part: frame	17 kg	20 kg
Mass heaviest component: leg support	2.4 kg	3 kg
Static stability downhill	-	13 °
Static stability uphill***	8°	15 °
Static stability sideways	15 °	15 °
Safe slope, use anti-tipper	-	7°
Seat plane angle	-3°	45 °
Effective seat depth	350 mm	500 mm
Effective seat width	350 mm	500 mm
Seat surface height at front with 7"	440 mm	440 mm
Backrest angle	90°	125°
Backrest height ex cushion from seat plate, max with back extender	500 mm	660 mm
Leg support foot plate to seat distance different adjustment profiles	261 mm	500 mm
Leg to seat surface angle	90 °	180 °
Arm support to seat distance	273 mm	329 mm
Front location of arm support structure - adjustable	300 mm	430 mm
Horizontal axle location	-10 mm	95 mm
Parking brake max slope	-	7°
Minimum turning radius, vertical leg supports	666 mm	702 mm

** For **Coyote T+** with seat width 500 mm, the overall width exceeds recommended 700 mm.

*** Least stable and most stable refers to the positioning of the anti-tippers. Always use anti-tippers driving uphill.

2. QUICK REFERENCE

-  The quick reference is not a replacement for the manual, only a reminder/ check list.
-  Adjust the wheelchair to the user: Adjust seat depth, foot plate height, armrest height, head support height and depth, chair back cushion height and eventually the wheelchair balance,

For more information about adapting the wheelchair to the user please see:
www.lohmedical.com

-  ANTI-TIP will secure the chair from tipping backwards as long as they are in use.
-  If in doubt - contact your dealer!
-  Drive carefully! Be sure to lock all handles properly.
-  Watch out for pinching danger when folding and unfolding, tilting, reclining and all other adjustment movements.
-  Never stand on the foot plates due to risk of tipping forwards.
-  The anti-tips are always used for the safety of the user.
-  When the chair is tilted rearwards the anti-tips must be activated. The brakes must be locked when leaving the user in the backwards tilted position.
-  Never lift the wheelchair by the leg supports, arm supports or head support.
-  Be aware that friction against push-rims can create a warm surface.
-  Surface temperature of metal parts in frame structure might increase when exposed to direct sunlight.
-  Salt water can increase risk of corrosion. Further precautions related to environmental conditions not needed.
-  If the chair has pneumatic tires: Make sure to check tyre pressure every week and inflate to keep 7" wheels at 36 PSI.

3. DESCRIPTION

-  **If any of these parts are missing, please contact your dealer.**
-  **For complete information, please contact your dealer.**
-  **Be aware that specifications may vary between countries.**
-  **Be aware that illustrations and pictures may differ from your received product.**

3.1 Information

Coyote T+ is an advanced mobility aid incorporating a range of adjustment solutions allowing optimum adjustment to the user:

Configurations - 45 degree tilt:

Front wheels	Rear wheels	seat height
7"	16"	435 mm

Dependent on wheel position the seat height varies.

Seat width:

Available by placing order:

350, 380, 400, 430, 450, 500 mm.

Seat depth:

Seat depth can be adjusted by repositioning of the chairback brackets measured from front seat plate to back rest cushions in steps: 350, 400, 425, 450, 475 & 500 mm. The main wheels are to be moved backwards in same steps to keep good balance in the chair.

Tilt range:

Coyote T+ wheelchair frame allows for extensive tilt range from 3° forwards to 45° (16" main wheels) backwards see configuration tabel to the left.

Recline range:

Coyote T+ wheelchair back allows for the back recline from 90° to 125°. When seat depth is 500 mm and the wheelchair has been tilted maximal 45°, there will only be 28° recline possible due to conflicts beteen arm support and wheels.

Manoeuvrability:

Coyote T+ has a wheel base where the foot print - and thereby the stability, balance and the ease of driving (manoeuvrability of the wheelchair) can be adjusted by adjusting the position of the rear wheel brackets, range -125 mm.

Coyote T+ is a relatively narrow chair - only 210 mm wider than the seat width, ensuring easy manoeuvring in narrow spaces.

4. FEATURES

4.1 Accessories


Not all accessories are available for all wheelchair configurations. Please check the home-page or ask customer service for further details.

STANDARD

Seat Cushion & Backrest Not included

Head support Comfort pressure distributive.


Pelvic guide - also to be used as knee support.

Seat plate aluminium

The seat can be pulled out and pushed in while the user sits in the seat. This eases transfer and positioning of the user. NOT to be used for seat depth extension.


ABDUCTION BLOCK

The block reduces abduction.

Small: 80 mm width

Medium 110 mm width

Large: 140 mm width

Universal leg supports with foot plates and heel strap height adjustable.


Height and depth adjustable arm supports Armpad 385 x 58 mm

REAR WHEELS

16" with drum brake Puncture proof Flexel or PU.

FRONT WHEELS / CASTORS

Available in 7"

TOOL KIT

MOUNTING OF HIP BELT


- Pull the belt thorough the hole in the hip belt bracket.

- Thread the belt back through the belt clamp.

- Fix the hip belt bracket to the rearmost hole in the back hinge, using the enclosed screws and nuts.

OPTIONAL (with extra cost)

5. ASSEMBLING AND ADJUSTING

5.1 Unpacking

See chapter 3 for details for required components (See chapter 5 & 6 for details for mounting).

1. Unpack all the parts, and check that everything is delivered according to the packing list.
2. Mount main wheels and front castors.
3. Mount the back rest gas cylinder, arm supports, cushions, head support and leg supports.
4. Mount accessories.

Weight of components (450 mm chair width):

Drive wheels: 1.7 kg each, Front castors: 0.9 kg each

Leg support Universal left & right: 1.3 kg (ea).

Necessary tools are described under each chapter.

Accessories described in chapter 4 are a presentation of options, and will be delivered with separate mounting descriptions.

When adjustments of seating and wheel positions are done in the defined positions by standard equipment, the adjustments will not exceed safe limits.

5.2 Initial adjustments to the user - overview

Start with user sitting in a normal upright position.

SEAT ADJUSTMENTS

It is vital to adjust the seat correct to give the user the support for his /her whole body.

The seat can be adjusted in depth and height.

The seat width is defined when ordering the wheelchair.


Adjustments to the seat may change the centre of gravity which effects the balance of the wheelchair.

Proper adjusted seat provide stability and the ability to manoeuvre the wheelchair safely.

A well balanced chair is easy to drive without easily tipping backwards.

Correct seat width depends on the user's hip width. It is measured while seating.

- The seat width is defining the wheelchair size when ordering.

Correct seat depth

depends on the user's thigh length. It is measured while seating - see the red horizontal arrow on the illustration.

ADJUSTING THE SEAT DEPTH

The seat depth can be adjusted both in the rear and in the front of the seat. Adjustments are done to give the user a comfortable seating position with proper lumbar support while the knee joint is aligned with the leg support knee joint. By adjusting the seat depth the chair balance and the driving characteristics may change.

Always start by adjusting the seat depth rearwards. Adjust the position of the recline gas cylinder to keep correct chairback angle. Adjust the main wheel position accordingly to seat depth. Secondly adjust the leg support fixing bracket to align the leg support knee joint with the user knee joint.

The seat depth can be adjusted 100 mm in the rear and 100 mm in front.

Correct seat height

The seat height is to be adjusted according to how the wheelchair shall be used:

- When driving with the push rims or being driven by an attendant the seat height is to be adjusted along with the leg supports. The leg support foot plates must have ca 30 mm free space underneath.

- When foot propelling - the seat height shall be adjusted for the user to have his/her feet flat on the floor.

Seat height depends on the user's calf length and is measured while seating - see the dotted vertical F-dimension on the illustration.

-  Please note: the leg support adjustments are described in detail in chapter 5.9.
-  Please, see the drawings on next pages regarding seat depth adjustments.
-  The relationship between seat angle and backrest angle should be kept constant when the tilt function is used to create variation for the user.

5.3 Wheels assembling

5.3.1 16" Drive wheels

To mount the 16" drive wheels the main wheel bracket must be mounted downwards:

- remove the small outer wheel cap,
- insert the bolt into the wheel bracket
- tighten the nut and bolt firmly.
- Repeat on the other side.
- Make sure the wheel has identical position on each side.

-  **2 pc 19 mm open-end spanner**

5.3.2 Drive wheels with quick release

-  To check that the drive wheel is properly attached to the hub, remove the finger from the central knob and pull the main wheel.
-  If the drive wheel does not lock, see chapter 5.3.3 and adjust.
-  If the drive wheel still does not lock, do not use the wheelchair but contact your dealer.
-  Sand and sea water (salt used for gritting in the winter) can damage the bearings of the main wheels. Clean the wheelchair thoroughly after exposure.

Choice of main wheel size and position in the wheel bracket - see table - allows for different driving characteristics such as hand propelling or attendant drive and seat height adjustments.

5.3.3 Correct frame angle

To ensure good driving characteristics you need to check and adjust - with drive wheels and front castors mounted - that the lower frame is horizontal and the front castor bearing house is vertical.

See the red lines.

Reposition the wheels until the frame is horizontal.

5.3.4 Main wheel quick release

The quick release axle links wheel and frame. The nut sitting under the head of the quick release axle can be adjusted to fit the axle tightly when inserted into the frame.

Adjust the nut backwards towards the axle head if the axle does not click in .

Adjust the nut forwards if the axle sits loose when inserted through wheel and brake.

2 pc 24 mm open-end spanner

5.3.5 Repositioning of the main wheel brackets for- and backwards

The wheel position (factory setting) is set at a point where the wheelchair has good driving characteristics and is stable even when tilting. It is possible to adjust the wheel distance by repositioning/moving the main wheel bracket backwards or forwards on the frame.

MAIN WHEEL BRACKET

- Remove the main wheel before adjusting.
- Loosen the 6 M8 screws holding the main wheel bracket (around the frame on both inside and outside).
- Slide the brackets with crossbar and anti-tippers (blue parts) forwards or backward to desired position.
- Tighten the screws well.

• Make sure the position is identical on the right and the left side of the wheel-frame; positions marked on the frame makes it easy to identify.

1 pc 6 mm Allen key

• The anti-tips follow the position of the main wheels and stays in correct position when the main wheel bracket is moved.


A scale on the wheel frame shows which wheel position is correct for the adjusted seat depth following the user needs.

Let the wheelchair roll forwards and check that it keeps a straight line. If it waver off to one direction, the main wheels position need adjustment.

The risk for tipping increases when the main wheels are moved forwards on the frame without also shortening the seat depth on the back brackets.


5.3.6 Height adjustments by Repositioning wheels in the main Wheels brackets

By repositioning the wheels in the main wheels brackets you can achieve higher or lower seat height - from lowest to heighest - by 2 x 25 mm.

2 pc 19 mm open-end spanner

Make sure to adjust the brake plate fixing as you reposition the wheels.


1 pc 5 mm Allen key

Remove the 4 hex button socket screws on each side holding the crossbar assembly to the wheelframe. See illustration.

1 pc 3 mm Allen key

Slide off the cross bar and anti-tipping device assembly from the wheel frame.

Loosen the wheel fixing brackets from the frame by unscrewing the 3 pc M8 hex socket screws - on each side, and slide them off the frame.

1 pc 6 mm Allen key

Slide the crossbar assembly with anti-tips on to the wheel frame.

Tighten the 3 pc M8 hex socket heads screws well on the wheel brackets - on each side.

Adjust the tramp pedals and fix them well to the frame with the 2 pc M5 hex socket button screws.

Check the position of the anti tips.

Check that the position of the wheel brackets are identical on both sides. A scale on the wheel frame shows which wheel position is correct for the adjusted seat depth following the user needs.

Adjust the user brakes according to the wheel position and size. Chapter 5.3.11 and 5.3.12

Let the wheelchair roll forwards and check that it keeps a straight line. If it waver off to one direction, the main wheels position need adjustment. The risk for tipping increases when

the main wheels are moved forwards on the frame without also shortening the seat depth on the back brackets.

5.3.7 Anti-tips

Coyote T+ has foot operated swing-away antitips. They follow the main wheel position and will normally be in correct position.

How to operate:

- Push the red anti-tip button downwards.
- Swing the anti-tip inwards, it will snap into position.

To swing it out

- Push down a little and pull it backwards. It will snap into active position.

Coyote T+ anti-tips can also be swung up - hand operated - when free space between main wheels is required.

How to operate:

- Pull the anti-tip backwards and swing up. It will snap into position.

ANTI-TIP HEIGHT ADJUSTMENT

The anti-tips can be height adjusted by repositioning the end piece. Remove the 2 screws and position in required holes. Height adjustment 25 mm

1 pc 4 mm Allen key

Anti-tip position depending on seat height:

Seat height 410 - anti-tip shortest position

Seat height 435 - anti-tip middle position

Seat height 460 - anti-tip longest position

- see illustration right.

Always check the anti-tip height. It may be necessary to adjust to get desired distance from anti-tip wheel to the ground - ca 30 mm.

5.3.8 Stepper pedal adjustment

Coyote T+ has a stepper pedal for the ease of tipping the wheelchair when passing obstacles.

The pedal is positioned according to the main wheel position. When the main wheel position has been adjusted it is recommended to adjust the length of the tramp pedal. It can be moved up to 100 mm.

- Loosen the M5 screw on the frame - see picture. Do not unscrew this screw!
- Position to desired length.
- Tighten the screw well.

1 pc 3 mm Allen key

Position the stepper pedal with equal length on both sides. A scale on the horizontal tube helps identify position.

5.3.9 Front fork

Front forks come as standard with quick release axles. The front fork is easily removed by pressing the quick release button on top of the bearing house.

5.3.10 Castor position in front fork

Mount the front castor in the hole of the front fork where the wheelchair frame is horizontal. Thereby the front fork is vertical and give good driving quality.

- Repeat on the other side.

2 pc 13 mm open-end spanner.

Choice of drive wheel size and position requires adjustments of the front castor position in the front fork.

5.3.11 Front casters

Removal

- Press the release button on top of the front fork bearing house - by pressing on the silicon cap covering the release button.

Mounting

- Lead the quick release axle into the bearing house and into the front fork.
- Pull the fork slightly to ensure that the fork is fully locked.

Sand and sea water (salt used for gritting in the winter) can damage the bearings of the front castors. Clean the wheelchair thoroughly after exposure.

5.3.12 The brakes


The **Coyote T+ 16"** wheels is fitted with

- attendant running drum- and parking brakes.

Operating and applying the attendant brake

The 2 hub brake handles for the drum brakes, which allow regulation of speed on hills and whilst travelling along are located on the push handles.

- To apply the brakes, pull the brake levers evenly and smoothly towards the handles.
- To activate the drum brakes for parking, push the little lever on the brake handle downwards while you hold the brake handle pulled.
- Next time you pull the brake lever the parking function is release.

The user parking brakes cannot be used as running brakes.

Do not leave the user in the wheelchair without activating the parking brakes.

5.3.13 User brake adjustment

The user brake bracket is fixed sitting in a profile that is fixed to the wheel-frame with 2 screws.

The brake bracket position can be adjusted forwards- backwards by loosening the 2 M6 screws and slide the bracket.

The adjustment is necessary when the main wheel position is changed.

IF THE BRAKE DOES NOT BRAKE PROPERLY:


Adjust the brake-brackets sitting on the wheelframe backwards or forwards till the distance between the braking bolt (B) and the tire is ca 3-5 mm when the brake is not active.

Test that the braking function is sufficient without being too hard to activate for the user.

5 mm Allen key (A)

The brake handle can be adjusted by loosening the screw (C) and position as required.


3 mm Allen key (C)

5.3.14 Drum brake adjustment

If the drum brake does not brake properly: Adjust the wire on one or both sides, adjust the foot screw 2-4 rounds out. Then re-check the brakes.

If the wire is too loose:

Adjust the foot screw all the way in. Tighten the wire by loosening the wire clamp before pulling the wire further through it. Tighten the wire clamp, and adjust the foot screw out again.

1 pc 10 mm open-end spanner.

To ensure the correct functions of the wire, these must never be taut.

5.4 Seat assembling - Wheelchair seat height

Correct seat height is based on a combination of identified intended use (foot propelling or hand propelling) and the size of the user (leg length).

5.4.1 Seat height

The seat height depends on:

- The size of the main driving wheels.
- The position of the driving wheels in the main wheel brackets.
- Size of front castors.
- Position of front castors and size of the front fork.

The Seat plate is an advanced seating platform developed to give persons sitting long time support and stability.

5.4.2. Seat depth adjustments

Seat depth adjustments are done to give the user a comfortable seating position with support under the thighs while having proper lumbar support and the knee joints are aligned with the leg support knee joint.

The seat depth can be adjusted both in the back and in the front of the seat. By adjusting the seat depth the chair balance and the driving characteristics may change. **A well balanced chair is easy to drive without easily tipping backwards.** Always start with adjusting the seat depth backwards, then adjust the leg support fixing bracket to make the leg support knee joint align with the user knee joint. If necessary the driving wheel position also needs to be change.

The seat depth can be adjusted 100 mm in the rear and 100 mm in front.

ADJUSTING SEAT DEPTH AT THE REAR BY MOVING THE CHAIRBACK HINGES.

The back rest hinge has 5 holes with 25 mm distance between them.

- The seat depth can be adjusted from 400–500 mm (Measured from front of seat plate to backrest tube without backrest cushion).
- On the frame under the chair plate there is mounted an extra locking head for the gas strut.
- The short locking head covers seat depth 400 mm, 425 mm and 450 mm.

- The long locking head covers seat depth 450 mm, 475 mm and 500 mm. (There is also a long back rest hinge which covers longer seat depths. See chapter 5).


Changing the seat plate size

- Loosen 4 pc M6 screws fixing the seat-plate part to the seat frame and replace with a size that fit the user.
- Tighten all screws well.

5 mm Allen key

By adjusting the seat depth, the chair balance and the driving characteristics may change. A well balanced chair is easy to drive without easily tipping backwards.

Chair width	Seat depth standar	Seat depth range
350 mm	380 mm	380 - 455
400 mm	405 mm	380 - 550
450 mm	480 mm	430 - 550
500 mm	480 mm	430 - 550

- When changing the position of the back rest hinge, also remember to change the position of the gas strut under the chair. The bracket has 3 holes. When the back rest hinge is set in the shortest position, the gas spring is placed in the front hole using the short locking head.

For each position the back rest hinge is pulled out, the gas spring is moved one hole towards the rear. When the gas spring reaches the rearmost hole with short locking head, change to long locking head and start from the front hole again.

When using the comfort seat plate it is recommended to change seat plate when changing seat depth. The seat plate comes in depth: 400 mm, 430 mm & 460 mm.

Seat depth in the rear is adjusted as follows:

- Release the tilt gas strut by tilting the seat all the way forward.
- Pull out the locking bolt for the backrest, and place the back rest forward in the seat.
- Remove the screws holding the back rest hinge, and set the hinge to a position which gives the backrest the required position. Do both sides at the same time.
- Place the screws back and fix them.

- Also remember to move the gas strut as described earlier in this chapter – to reach 90° back angle when the gas strut bottoms.

When adjusting the seat depth at the rear, the screw in the backrest hinge must be tightened with 14 Nm.



**6 mm Allen key.
13 mm open-end spanner.**

If the user requires another back rest angle than what is standard, it is possible to change the position of the gas strut in three positions under the chair in front.

When changing seat depth, you also change the tipping point of the chair. This can be prevented by changing the position of the main wheel brackets.

ADJUSTING SEAT DEPTH AT THE FRONT

It is possible to adjust the seat depth up to 90 mm at the front. For easy seat depth adjustment seat depth extension brackets can be delivered (+45 mm). They sit on the wheel frame. By restless users the extension piece should not be pulled out more than 50 mm.

Do the following:

- Loosen the screw holding the extension piece for the leg support.
- Set the extension piece to the required position.
- Place the seat plate extender on the leg support extension pieces.
- Fix the screws, tighten them with 25 Nm.

6 mm Allen key.

By setting the leg support extension pieces to different positions – different position / depth of the leg supports, it is possible to compensate for a rotated pelvis or different length of thighs.


Never use the sliding seat plate in forward position as seat depth extension.

5.5 Chair back assembling and positioning

To mount the recline gas spring:

- Lift the back rest by its push handle with one hand,
- Lead the gas spring locking head into the plastic bracket which is sitting on the lower cross beam at the rear.
- Secure the back rest by pushing the locking bolt in from the side, through the bracket and gas spring locking head.

To check that the back rest is locked, grip the push handle and press the back rest forward.

If the backrest falls forward - repeat the locking procedure or contact your dealer.

5.6. Assessment and choice of back support gas spring

The chair back position / recline angle is defined by the chair back gas spring. There are different spring strengths available for the **Coyote T+**:

INDICATIVE correlation between body weight and suitable recline gas spring force: Standard delivered is 250 N recline gas cylinder. Please contact your dealer if you need other strength.

5.6.1 Change of recline gas spring

- Remove back cushion.
- Remove the arm supports.
- Remove the back spring splint.
- Fold the chair-back forwards.
- Unhook the wires attached to the cylinder head.
- Remove the gas spring by unscrewing it from the cylinder head.
- Replace the gas spring with a new - with correct strength for the user.
- Mount all parts that was removed.
- Adjust the wires.
- Fold up the chair-back and test if the solution suits the user strength.

5.6.2 Change of tilt gas spring

The **Coyote T+** is equipped with 2 tilt gas springs. Standard delivered tilt gas spring strength is 300N. Replacement of tilt gas spring:

- Unscrew the nut on the upper end of the spring.
- Loosen the nut at the lower end of the spring.
- Unscrew the gas spring.
- Replace with the new spring.

Screw the new spring into the spring head and then into the spring bracket.

- Tighten the nuts on the bottom and on the top.
- Check that the tilt function is working correct.

1 pc open end spanner

5.7 Head support assembling and positioning

5.7.1 Head support

A - Lever for depth adjustment

B - Wheel for angle adjustment

C - Lever for height adjustment


D - Head support bracket.

- Place the squared nut in the trace of the head support bracket as shown below.

- Place the head support in the head support bracket.

- The height and the depth of the head support is set to the required positions and tightened.

- The head support bracket is fixed by tightening the four screws two by two diagonally so the bracket is fixed with the same strength divided on the four screws.

MOUNTING AND REMOVING

- To mount the head support, loosen the red italian grip on the push bar.

- Insert the head support rod into the head support bracket.

- Fix at any desired height by tightening the italian grip

HEIGHT ADJUSTMENT OF THE HEAD SUPPORT

The middle of the head support cushion shall be positioned to meet the head slightly lower than the most rearward part of the head.

DEPTH ADJUSTMENT OF THE HEAD SUPPORT

To be adjusted by the therapist. Loosening the Italian grip on the side allows for the head support cushion to move backwards and forwards.

ANGLE ADJUSTMENT OF THE HEAD SUPPORT

Loosening the star-wheel at the end of the horizontal bar allows for adjustment of the angle of the head support cushion.

Adjusting the head support sideways:

- The head support adapter on the push bar can be moved both to the right and left, giving the possibility to accommodate special needs for head support.
- Loosen the four screws holding the adapter together.
- Move the adapter to the required position and fix the adapter by tightening the screws diagonally.

Remember to release the levers when adjusting the head support.

If the head support stand does not fit the bracket perfectly, the bracket is probably fixed too tight or unevenly. After fitting the head support, fix it properly by tightening the little set screw in the centre on top of the head support bracket using an Allen-key.

If the head support seems too short in height, it can be turned 180° by releasing the adjustment wheel at the rear of the horizontal bar (B).

5.8 Arm supports adjusting

5.8.1 Arm supports height

Adjust the arm support height to fit the user. The height should give good underarm support without lifting the shoulder. There are 5 height positions available, with 15 mm distance between each.

- Unscrew and position the arm support in correct height for the user.

5.8.2 Arm supports depth

If the adjustments at the front of the arm support is not giving you enough support the arm pad can be moved forwards.

- Loosen and slide the arm pad to desired position.

1 pc screwdriver for M5 mm Phillips screw

ARM SUPPORTS PRECISION -

If the arm support vertical profile fits either too loose or too tight in the arm support bracket it can be adjusted to better fit.

- Tighten or loosen the 2 M8 set screws sitting on the rear side of the arm support brackets.
- Test the fit by inserting the arm support.

5.8.3 Arm supports handling

The arm support is inserted into the arm support bracket. It is removed by lifting it up and out of the bracket.

When the arm support is removed there is free space for side transfer.

With the arm supports removed there is good access to place sails etc behind the back of the user.

Always lock the brakes when you plan for actions and remove the arm support.

5.9 Leg supports

Coyote T+ is as a standard delivered with Universal leg supports.

Coyote T+ allows for adjusting the leg support forwards. Extra cost, request them in order.

5.9.1 Leg support height and length adjustments

Leg support adjustment takes place on 2 different locations.

First adjustment is the foot plate height. The second adjusting the seat plate depth and the seating unit position in the frame. The goal is having the leg support knee joint and the user knee joint aligned in the same axis.

Measure the user's lower leg length F - 90 degrees bent knee - measure from underside thigh to the underside heel of normal used shoe. The knee joint of the leg support may be hidden under the knee upholstery. Adjust the foot board height according to the lower leg length measured.

Make sure there is enough free space under the foot-plates for the wheelchair to pass minor obstacles.

Slightly tilting the seating unit may help.

If not sufficient, the chair height must be adjusted as described in this user manual chapter 5.3

CORRECT SEAT DEPTH:

Adjust the seat plate and position as described in chapter 5.4.

- The foot plates are hinged, and can be angled in fixed positions.

5.9.2 Universal leg supports

The Universal leg supports are swing-able, height adjustable and removable. They are adjusted to a fixed angle. Heel straps prevent the feet from sliding from the foot plates.

UNIVERSAL LEG SUPPORT MOUNTING INSTRUCTION:

- Fold the foot plates up.
- Hold the leg support on the top joint, and place it in the leg support bracket on the wheelchair in an outward turned angle as shown in the picture.
- Swing the leg support inwards and push slightly downwards until it goes into locked position.


ANGLE ADJUSTMENT

- Loosen the star wheel on the outside of the leg support.
- Lift the lower part of the leg support to required angle.
- Tighten the star wheel.

FOOT PLATE HEIGHT ADJUSTMENT:

The foot plates are step less height adjustable.

- Unfix the adjustment screw so that the adjustment bar moves freely.
- Slide the foot plate to required height,
- Tighten the screw.

5 mm Allen key

FOOT PLATE ANGLE ADJUSTMENT:

- Loosen the screw on the outside of the foot plate - see picture - with an Allen key.
- Tilt the foot plate to required angle and tighten the screw.

6 mm Allen key

Never stand on the foot plates due to the risk of tipping forward.

5.9.3 Angle adjustable leg supports

The angle adjustable leg supports are swing-able, height adjustable and removable. They come with height- and depth adjustable calf supports. The foot plates are hinged, and can be angled in fixed positions.

The foot plates come with a lock connecting the 2 plates which makes the plates stronger. If locking is not wanted, the bolt can be removed by using an Allen key with no loss of functionality.

ANGLE ADJUSTABLE LEG SUPPORT MOUNTING INSTRUCTION:

- Fold the foot plates up.
- Hold the leg support on the top joint, and place it in the leg support bracket on the wheelchair in an outward turned angle as shown in the picture.
- Swing the leg support inwards and push slightly downwards until it goes into locked position.

ANGLE ADJUSTMENT

- Loosen the star wheel on the outside of the leg support.
- Lift the lower part of the leg support to required angle.
- Tighten the star wheel.

Be aware of squeeze hazard between moving parts.

FOOT PLATE HEIGHT ADJUSTMENT:

The foot plates are step less height adjustable.

- Unfix the adjustment screw so that the adjustment bar moves freely.
- Slide the foot plate to required height,
- Tighten the screw.


As an accessory a star knob or a lever handle can replace the screw.

FOOT PLATE ANGLE ADJUSTMENT:

- Loosen the screw on the outside of the foot plate - see picture - with an Allen key.
- Tilt the foot plate to required angle and tighten the screw.

5 mm Allen key

LOCKING AND RELEASING THE FOOT PLATES

- To lock the foot plates, let the right foot plate fall over the bolt standing out from the left one, it clicks into lock.

- To release the foot plate pull the plastic lock between the foot plates and lift the right foot plate up.

While making the adjustment, there must be no load on the foot plates.

For outdoor use, there should be a clearance of 40-50 mm between the foot plate and the ground.

Never stand on the foot plates due to the risk of tipping forwards.

LOCKING AND RELEASING THE FOOT PLATES

- The foot plates come with a lock connecting the 2 foot plates which makes the plates stronger.

- To lock the foot plates let the right foot plate fall cover the bolt standing out from the left one, it clicks into lock.

- To release the foot plate push the plastic knob under the right foot plate and lift the right foot plate up.


While making the adjustment, there must be no load on the foot plates.

For outdoor use, there should be a clearance of 40-50 mm between the foot plate and the ground.

When adjusting leg support angle, be aware of squeeze hazard between moving parts.

5.9.4 Leg support mounting and removal

Mount the leg supports to the wheelchair. Ensure an open angle of the leg support when inserting:

Position the bolt vertical into the frame hole turning it ca 30 degrees outward for it to enter easy. Turn it inward until it clicks into user position.

When mounting the leg supports, always check that there are no conflict between the front castors and the foot plates. Adjust the leg support angle if necessary.

5.9.5 Foot plate lock

Leg supports can have foot-plates replaced by a foot-board. With both leg supports in place, the foot-board folds down into the lock on the right leg support.

The lock is closed by turning the outer ring 90 degrees forward!

The foot-board lock **MUST** always be closed when the chair is occupied.

5.10 Abduction block

Coyote T+ seat plate is prepared for mounting abduction block to the front of it.

- Screw the abduction block bracket to the seat plate in the front.
- Position it in the centre of the plate.

Adjust the height to fit the block tightly over the seat cushion.

4 mm Allen key

6. EVERYDAY USE

Coyote T+ wheelchair is an advanced mobility aid. Please use some time to study this manual and become familiar with the wheelchair and all its possibilities.

Max user weight is 136 kg.

After the wheelchair has been correctly adjusted to the user described in the previous chapters, you as an attendant need to learn how take advance of the tilt and recline features to give the user the most benefits of the chair.

SEAT UNIT ANGLE - TILT

The seat unit angle is regulated with the tilt handle on the left push handle. The seat unit can be tilted from 0° to +45° when 16" wheels are mounted.

The tilt function is used for giving variation of the seating position for the user.

BACKREST ANGLE - RECLINE

The backrest angle is regulated with the recline handle on the right push handle. The recline angle can be regulated from 87° - 125°.

6.1 Key words regarding tilt and recline of comfort wheelchairs

Tilt and recline are the basic benefits of a comfort wheelchair. It allows for varying seating positions during the time in the wheelchair, for creating correct position, maintaining stability, stimulate ADL and protecting skin.


We have reviewed the clinical evidences regarding tilt and recline, and found there are several studies and best practice guidelines suggesting that the tilt and recline sequence is important to reduce shear and sliding:

When changing to a resting backwards leaning position: tilt first, then eventual recline.

When bringing the user upright from a resting position,: the sequence should be reverse; first correct the recline angle then the tilt. The most shear is induced when going upright from a recline and tilted position.

The recline angle is to be adjusted by the therapist to accommodate the user's best seating position, Always returned to, this angle when sitting in upright seating position.

6.2 Decrease the possibility of sliding, shear and pressure sores

By upright normal seating the muscle tone of the neck and back should be as low as possible for the user to prevent sliding. The therapist will identify the correct recline angle. A change of the recline angle from this position can disturb the correct body position, and cause an increased muscle tone in the neck.


If the recline function is used during a transfer situation or other situations, it is very important that the recline angle is adjusted back to the correct, original position when the user is back to a normal seating position.

Wrong use of recline can cause loss of position, an increased risk of sliding, an increased danger of shear strain and pressure sores.

6.3 Operating recline: reclining the chair back

Pull the recline button on the right side of the push handle and recline the chair back. It will stop where you let loose.


MAKE SURE THAT THE USER IS SAFE WHEN TILT OR RECLINE ARE GOING TO BE ADJUSTED:

The tilt and recline functions are «one hand operations». This is a benefit for the user: The care giver is able to establish eye contact and communicate with the user when the tilt or recline function is going to be used. The user can feel more safe when tilted or reclined while eye contact is maintained between attendant and user.

6.4 Operating tilt: tilting the seating unit


Pull the tilt handle on the left side of the push handle and tilt the seating unit to desired position. The relative angle between the body upper and lower parts remain the same when the seating unit is tilted. **Coyote T+** can be forward tilted -3° -and backwards tilted $+30^{\circ}$ (45° with 16" main wheels). A rearwards tilted seat unit can create a position, with decreased risk of sliding, an decreased danger of shear strain and pressure sor

A forward tilted seat unit brings the user in a position where activities, for instance by a table or by standing up from the wheelchair, are supported.

Never leave the user alone when the seat is tilted forwards. The user may slide forwards out of the chair. 0° - 30°

Make sure the belts and harnesses are fixed in the best way for the user.

A hip belt is of benefit for the user to keep the seating position and for the safety.

Upper body harnesses can be useful to connect the user to the back support. When using upper body harnesses the hip belt must always be locked first.

Ankle huggers can be beneficial by connecting the feet to the foot plate- helping preventing them from sliding off.

6.5 Practice the use of the Coyote T+ tilt and recline system together with the user

- Practice using the tilt function leaning the complete seating unit backwards and forward, experiencing and testing how the user react.
- Tilt is recommended to be used for varying seating positions.
- Practice using the recline function by changing the back angle only. This is helpful when using a hoist to transfer the user in and out of the chair. Please return to correct recline angle after transfer.

6.6 Transfer to and from the wheelchair

Techniques for transferring to / from the wheelchair should be practiced well with the persons involved.

Some important advices for preparation of the chair:

- The wheelchair should be placed as close as possible to the destination of the transfer.
- Make sure the brakes of the wheelchair are activated to prevent it rolling away.
- Remove the leg support.
- Lift away the arm support on the side of the transfer.
- Eventually remove the head support.

6.7 Using a patient hoist

- Make sure the brakes of the wheelchair are activated to prevent it rolling away.
- Tilt the chair slightly backwards.
- Remove the leg support.

- Open the back support angle slightly while holding the user and adjusting the hoist behind the users back.
- Option: Remove the arm supports to get closer to the patient and / or remove the head support.
- Replace the components when the transfer is finished.

6.8 Daily attendant handling

6.8.1 Push bow

The push bow are easily height adjusted at any time enabling good control over the wheelchair in the varying positions.

For safe manoeuvring the wheelchair, and preventing injuries for the attendant, the push bow should be adjusted according to attendants arm height.

6.8.2 Anti-tips

Coyote T+ has foot operated swing-away anti-tips to ensure a safe and stable wheelchair. They are premounted according to standard set up. Press the button down while you swing the antitip inwards or outwards.

The anti-tips can also be swung up to give free space between the wheels. Pull backwards and swing up.

6.8.3 Tramping pedal

The tramping pedal is used for tipping the wheelchair backwards when passing obstacles.

- Step with one foot on the tramping pedal.
- Pull and slightly press the push handles downwards and at the same time.
- Tip the wheelchair.

6.8.4 Brakes

Coyote T+ with 16" main wheels is equipped with drum brakes that are operated with the brake handles on the push handles.

The brake handles on the push handles can be used as running brakes like on a bicycle.

ATTENDANT BRAKES HAVE PARKING BRAKES

- To apply the brakes, pull the brake levers (1) evenly and smoothly towards the handles.
- To activate the drum brakes for parking, push the little lever (2) on the brake handle downwards while you hold the brake handle pulled.
- Next time you pull the brake lever the parking function is release.

Always use the parking brake when the wheelchair is left without attendant standing close. Lock the parking brakes on both sides.

6.8.5 Arm supports

Coyote T+ is equipped with removable arm supports.

LOCKING - UNLOCKING

The arm supports are simply lifted out of the arm support bracket.

The high position screw keep the arm support at correct height every time.

Mounting the arm support:
inserting the pole into the arm support bracket.

6.8.6 Head support-mounting

Coyote T+ is a standard fitted with head support.

Head support adjustments:

The head support adjustment should be set by the therapist adjusting the wheelchair to the user.

The height of the head support needs to be adjusted to fit directly behind the head giving support to the lower part of the scull.

The depth of the head support should be adjusted to barely touch the back of the user's head when sitting relaxed.

HEAD SUPPORT MOUNTING AND REMOVAL


The head support for is easily mounted and removed:

Open the italian grip 1. The head supports is released and can be lifted up and removed.

To mount - insert the vertical head support pole into the head support bracket and tighten well.

Coyote T+ can be equipped with Dynamic head support.

The dynamic head support follows the backwards movement of the neck. This is beneficial both for users with involuntary movements and users nodding or banging the head.

6.8.7 Leg supports - mounting

Leg supports removal:

- Unlock the foot plates - if a lock is mounted.

Pull the leg support slightly upwards.

Swing it outwards before lifting and detaching it from the chair.

Leg supports mounting:

Push the leg support fixing bolt vertical into the leg support fixing bracket on the wheelchair - keeping the leg support slightly turned outwards.

The leg support will fall in place and turn into place itself.

FOOT PLATE LOCK

Many of leg supports have foot plates with lock. The lock connects the foot-plates and help stabilize the foot-plates.

If a footplate lock is mounted it must be opened before the foot-plates can fold up.

To open: pull the red button sitting in the gap between the foot plates and fold up.

When the foot-plates fall down they will most likely fall into the lock by themselves. Please check.

Always lock the foot board when the user sits in the chair. Failing to lock may lead to damage to the leg supports and eventual also harm to the user.

7. TRANSPORT

7.1 Transport in car

Whenever possible, transfer to a car seat with vehicle safety belts when you are travelling with a car. Secure the wheelchair or store it in the cargo area of the car.

Coyote T+ has been successfully crash tested in a forward facing orientation with both pelvic and shoulder belts, according to the requirements of ISO 7176-19:2022 and is approved to be used as a seat in a vehicle.

Coyote T+ is tested with a combined wheelchair and occupant restraint system W120/DISR from Unwin Safety Systems. Always use approved wheelchair and occupant restraint system (ISO 10542:2012) for fixing the wheelchair in the vehicle. Use a 4-point strap-type tiedowns to secure the wheelchair in the vehicle. The rating for the wheelchair's accommodation of vehicle anchored belt restraints is rated A = good.

REMOVE ACCESSORIES

Before using the **Coyote T+** as a seat in a car, be sure to remove and secure all auxiliary parts and accessories (e.g. trays and abduction block) that may fall off the chair in case of an accident and secure them safely elsewhere.

Coyote T+ has been crash tested without any power assistant device etc. If, at a later point of time a power kit, stair climber etc. is mounted, you need to check if your power assistant device is crash tested and approved for wheelchairs being used as seat in a car. If not, the assistant device must be dismantled and secured elsewhere, when the wheelchair is used as a seat in a car.

Coyote T+ seat width 350-500 mm is approved for user weight up to 136 kg when used as a seat in a vehicle.

SECURING THE WHEELCHAIR

Raise the chair to an upright seat position with max 10 degree tilt and 10 degree recline.

Swing the push bow up and fix it behind the head support.

For fixing the wheelchair in the vehicle only the designated securement points should be used.

The chair is marked with stickers showing the wheelchair securement points. Tiedown ear dimension 25x60mm.

In front: Use hook or strap attachment.

In the rear: Hook on strap or a hook/carbine hook into the car fixing loops on the frame.

The angle of the straps should be close to 45°

SECURING THE USER


Always use 3-point car occupant restraint belts.

Always use both pelvic and shoulder restraints to reduce the possibility of head and chest impact with vehicle components. Watch out that the belt is not twisted and the release buckle will not get in contact with the chair in case of a crash.

Make sure the car safety pelvic-belt lay tightly across or in front of the pelvis - the angle between pelvic belt and the horizontal between 30 - 75 degree, the steeper angle the better. The shoulder belt must lay close to the body of the user and not across arm supports, wheels etc. See illustration. The shoulder belt must not lay across arm supports, wheels etc. The corrective harnesses used in the wheelchair are not safety belts.

 **If the user is 1.85 m or longer, the back extension kit and a 500 mm vertical bar for the head support must be mounted when is used as a seat in a car.**

 **If a head support is mounted correctly it is very stable but does not replace the need for an external neck support mounted in the car.**

 Never use the wheelchair as a seat in a car if it has been involved in an accident with impact, before the chair has been inspected and approved for this by the manufacturer's representative.

 Alterations or substitutions should not be made to the wheelchair securement points or to structural and frame parts or components without consulting the wheelchair manufacturer.

 The weight of wheelchair cushions varies depending on the size and configuration of the wheelchair.

7.2 Folding for transport

When the wheelchair is unoccupied, fold it as described below. Put the wheelchair in the trunk or the back seat. When placed in the back seat, secure the frame using safety belts.

- Remove the head support (chap. 6.8.5)
- Swing push-handles inwards (chap. 6.8.1)
- Revolve the arm supports up (chap. 6.8.4)
- Remove leg supports (chap. 6.8.6)
- Release backrest and fold it forwards (chap. 5.5)
- Remove main wheel (chap. 5.3)
- Remove front castor (chap. 5.3.3)

7.3 Transport in air-plane

Coyote T+ may be transported in airplane without any restrictions.

Coyote T+ is equipped with 3 gas springs. These are however not classified as dangerous goods.

Coyote T+ is configured with tilt and recline functions which are activated with gas spring(s).

Contrary to general dangerous goods instruction UN3164, the IATA-DGR (special regulation A114) rules that the goods that contain gas and are determined to function as shock absorbers (including energy-absorbing devices or pneumatic springs) are NOT subject to the transport instructions i.e. they are indemnified from the following requirements:

- a) Each article has a gas volume which does not exceed 1,6 l and a charge pressure not exceeding 250 bar, where the product of the capacity expressed in litres and charge pressure expressed in bars does not exceed 80.
- b) Each article has a minimum burst pressure of 4 times the charge pressure at +20 degree Celsius for products not exceeding 0,5 l gas space capacity.
- c) Each article is made of material that will not fragment.
- d) Each article was manufactured in accordance to quality standard which is approved by the responsible national authority.
- e) It is proven and shown that the article relieves its pressure by means of a fire degradable seal or other pressure relief device such that the article will not fragment and the article does not rocket.

7.4 Travelling on public transport

The wheelchair should be put in a special area for wheelchairs.

The wheelchair should face opposite the direction of travel. The back of the wheelchair must be located against a fixed object such as a row of seats or a partition. Make sure the user can easily reach any hand rails or handles. Use belts and harnesses in the chair to hold the user.

Use safety belts if available to secure the user in the vehicle.

8. MANOEUVRING

8.1 General techniques sitting in your chair

Once on your chair, make sure that you are comfortably positioned and that the leg supports and arm supports have been adjusted to suit your needs and that you have good back support.

GUIDELINES FOR USE

Always drive carefully. Make sure your seating position is tilted and reclined in a way that you do have the overview when driving.

BALANCE

The weight load and balance of the chair influences the manoeuvring ability of the wheelchair. The weight, size and sitting position of the user are also influencing factors.

The size and position of the wheels will influence the driving performance.

The more weight placed over the main wheels, the easier it is to manoeuvre. If heavy weight is placed over the front castors, the chair will be heavy to manoeuvre.

-  Avoid carrying bags on the push handles, it may make the chair unstable.
-  The large available tilt and recline positions of the chair influences the manoeuvrability of the wheelchair. When large tilted/reclined position is used, the chair is not meant to be driven long distances. The attendant can drive the wheelchair with large tilted/reclined position for shorter distances like transferring from room to room.
-  **RIGIDITY**
The large available tilt and recline movements have influence on the rigidity of the wheelchair.
-  **PASSING OBSTACLES**
The anti-tippers may get in conflict with steps and obstacles when passing. Swing the antitippers away when passing an obstacle.
-  **STEP APPROACH:**
Always approach a step in slow motion preventing the front castors to hit the step with force. The user could fall out of the chair by the impact. The front castors could brake.
-  **DRIVING FORWARD DOWN STEPS/ SIDEWALKS**
Be cautious that you do not drive down steps higher than 30 mm. The leg supports may hit the ground first. Thereby you might loose the control and the leg supports may brake.
-  **PARKING:**
Increase the underneath support of the wheelchair by moving the chair about 100 mm backwards making the front castors turn forward.
-  **DRIVING ON SOFT, ROUGH OR SLIPPERY GROUND** can make safe manoeuvring more difficult as the wheels may loose traction and it is difficult to control the wheelchair.
-  Leaving the user unattended: If the user is left alone in the wheelchair, always lock the brakes and secure that the anti-tips are turned down.
-  Emergency escape routes:

Please be aware that wider wheelchairs have bigger turning radius and reduced manoeuvrability in vehicles. Smaller wheelchairs generally provide greater ease of vehicle access and manoeuvrability to a reach forward facing position.

8.2 Driving techniques - step up

The anti-tippers can be folded away when passing steps and obstacles. Be careful when you have to pass heights.

Companions, drive up a step forwards:

- Angle the wheelchair backwards on the drive wheels by stepping on the tipper pedal before pulling on the push-handles.
- Push forward till the front castors stand on the step and continue to push while you also lift holding on to the push-handles.

Users, drive up a step backwards:

This technique is only useful if the step is low. It also depends on the clearance between the foot plates and the ground.

- Drive the chair backwards towards the step.
- Make a firm grip on the push rims and move the body forward while pulling.

Companions, drive up a step backwards:

- Pull the chair backwards next to the step.
- Angle the wheelchair backwards, moving the front castors slightly up in the air.
- Pull the wheelchair up by the push handles onto the step and go backwards long enough to put down the front castors on the step.

8.3 Driving techniques - step down

The anti-tippers can be folded away when passing steps and obstacles. Be very careful when you have to pass small heights.

Companions, drive down a step forwards:

- Fold the push-bow down.
- Angle the wheelchair backwards, moving the front castors slightly up in the air.
- Drive carefully down the step and angle the wheelchair forward putting the front castors back on the ground.

Companions, drive down a step backwards:

- Move the wheelchair backwards to the step.
- Drive carefully down the step and move the wheelchair backwards on the main wheel until the front castors have come away from the step.
- Put the front castors down on the ground.

8.4 Driving techniques - slope

Important advise for driving down and up hill avoiding the risk of tipping.

Avoid turning the wheelchair in the middle of a slope. Always drive as straight as possible.

It is better to ask for assistance than taking risks.

Driving uphill:

Move the upper part of the body forwards in order to maintain the balance of the chair.

Driving downhill:

Move the upper part of the body backwards to maintain balance of the chair.

Control the speed of the chair by clutching the push rims. Do not use the brakes.

8.5 Driving techniques - up and down stairs

Due to the weight load - with wheelchair and user, we do not recommend driving a wheelchair up or down stairs.

The strain on assistants will be more than the manual handling regulation and the ergonomic assessment max load.

The risk for losing control during stair passing can lead to danger for both wheelchair user and assistants.

Use ramps and elevators when moving the user to higher or lower floor.

Never use escalators, even if assisted by a companion.

Do not lift the wheelchair holding onto the push-bow, leg-, arm- or head supports.

8.6 Transfer

Techniques for transferring to/from the wheelchair should be practiced well with the persons involved. Here, we give some important advices for preparation of the chair.

With or without companion - sideways.**Before transfer:**

- The wheelchair should be placed as close as possible to the destination of the transfer.
- Pull the wheelchair backwards 50-100 mm in order to make the front castors turn forward.
- Lock the brakes.
- Remove leg support and swing up the arm support on the side of the transfer.

With or without companion – forwards.**Before transfer:**

- The wheelchair should be placed as close as possible to the destination of the transfer.
- Pull the wheelchair backwards 50-100 mm in order to make the front castors turn forward.
- Lock the brakes.
- Tilt chair forward.

Using a lift / hoist:**Before transfer to chair:**

- Tilt the chair slightly backwards
- Remove the head support
- Remove the leg supports
- Open the back rest angle slightly
- Replace the components when transfer is finished.

Never stand on the foot plates due to the risk of tipping the chair forwards.

8.7 Point of balance

Adjust the point of balance by adjusting to correct seat depth for the user. This will move his/ hers centre of gravity and balance the wheelchair.

In addition change the position of the main wheel in the main wheel bracket.

- Move the main wheel hub and the main wheel.

(Chap. 5.3.1)

- Adjust the drum brakes. (Chap. 5.3.6)

- Check that the main wheel and quick-release are locked properly. (Chap 5.3)

2 pc 24 mm open-end spanner | 5 mm Allen key

The Wheel bracket allows for changing wheel position and thereby changing the seat height.

The point of balance is also influenced by adjusting the seat.

8.8 Lifting the wheelchair

The unfolded wheelchair without user should be lifted from 2 persons holding on to the frame and push bar only.

It is marked with a symbol where it should be lifted.

Never lift the wheelchair when holding on to the leg-, arm- or head supports. They may detach and the wheelchair may fall and get damaged.

Never lift the wheelchair with a user in it.

9. MAINTENANCE

9.1 Maintenance instructions

You as a user of the wheelchair (and your attendants and family) are responsible for the everyday maintenance of the chair. Clean it regularly. Do the maintenance to ensure safe and long time reliable functions and hygienic appearance.

Frequency	Weekly	Monthly
Check defects / damages e.g. breakage / missing parts	X	X
Washing of wheelchair		X
Washing of cushions		X
Check anti tip function		X
Check brake adjustment		X
Check tyre wear		X
Oiling of bearings with bicycle oil		X
Grease vertical leg support profiles with white Vaseline	X	

9.2 Cleaning and washing

1. Remove cushions before washing the wheelchair.
2. Clean the frame using water and a rag.
3. We recommend using soft soap.
4. Rinse the wheelchair well using clean water to remove all the soap.
5. Use methylated spirit to remove any dirt left.
6. Clean cushions and covers according to instructions printed on cushions.

DISINFECTION OF THE WHEELCHAIR

Remove cushions. See separate washing instruction above.

Wipe disinfection: use a soft rag wetted with Hydrogen peroxide or technical alcohol (isopropanol) and wipe the whole chair clean.

Check / re-adjust screws and nuts at regular intervals.

Sand and sea water (salt used for gritting in the winter) can damage the bearings of the front castors and main wheels. Clean the wheelchair thoroughly after use.

* As a rule of thumb, use oil on movable parts and all bearings. We recommends use of ordinary bicycle oil.

9.3 Long term storing

If the wheelchair is to be stored for longer time (longer than 4 months) no special actions are needed. We recommend that the chair is cleaned before storing. Before it being used again, complete the above maintenance instructions.

SPARE PARTS

Coyote T+ is built of modules. OneLoh carries stock of all parts and is ready to supply. Necessary instructions for mounting will follow the parts.

These parts can, if needed, also be removed and sent to manufacturer / distributor upon request.

Parts related to wheelchair frame construction must be handled by manufacturer or authorized service facility.

If defects or damages occur, please contact your dealer.

Original paint for repair of repair of scratched, can be ordered from Loh Medical.

10. TROUBLESHOOTING

Symptom	Reason /Action	Reference in Manual
The wheelchair is going askew	- The angle of the bearing house might not be 90°. - Check that the front castors are fitted in the same height. - The main wheel hubs might be incorrectly mounted. - One of the brakes might be too tight. - The user are sitting very askew in the chair. - The user might be stronger on one side than the other.	5.3.2 5.3.3 5.3.4 5.3.5
The wheelchair is heavy to manoeuvre	- The main wheel hubs might be incorrectly mounted. - Clean the front castors and forks for dirt. - Too much weight over the front castors. (Adjust the point of balance by moving the seating unit backwards or move the main wheel brackets more forwards).	5.3.2 5.3.3
The wheelchair is hard to turn	- Control that the front castors are not fixed too tight. - Too much weight over the front castors. (Adjust the point of balance by moving the seating unit backwards or move the main wheel brackets more forwards).	5.3.3 5.3.4
The front castors are wobbling	- The front castors are not fixed properly. - Check that the front forks are fitted in the same height. - The angle of the bearing house might not be 90°. - Too much weight over the front castors. (Adjust the point of balance by moving the seating unit backwards or move the main wheel brackets more forwards).	5.3.2 5.3.3 5.3.4
The main wheels are difficult to take off and put on	- Clean and grease the quick release. - Adjust the length of the hub bushing.	5.3.6
The brakes are not functioning well	Adjust the drum brake.	5.3.7
The wheelchair feels "shaky"	Check screws and adjustment points in general.	

 **Please contact your dealer for information about authorized service facilities that can give support if solution is not reached in this form.**

 **When in need of spare parts, please contact your dealer.**

 When making changes affecting frame construction, contact dealer / manufacturer for confirmation.

11. TEST & WARRANTY

11.1 Tests

Coyote T+ is tested and have been approved for both indoors and outdoors use.

The chair is CE marked.

MAXIMUM USER WEIGHT: 136 kg

**It is tested by
TÜV SÜD Product Service GmbH
according to DIN EN 12183:2014.**

Coyote T+ is crash tested at TASS Netherlands according to ISO 7176-19: 2022.

It is approved for being used as a seat in a vehicle.

Maximum user weight when used as seating a car: 136 kg

11.2 Warranty

Loh Medical is providing you with a 5-year warranty on all frame components and on the cross-tube assembly. There is a 2-year warranty on all other CE labelled components.

Loh Medical is not responsible for any damage resulting from inappropriate or unprofessional installation and / or repairs, neglect, wear, from changes in wheelchair assemblies or instructions not approved by Loh Medical or by use of spare parts delivered or produced by third parties. In such cases, this warranty shall be considered null and void.

This warranty is only valid when the user use, maintain and handle the wheelchairs as described in this user manual.

11.3 Claim

If a product has developed a fault during the warranty period as result of a defect in design or manufacturing, you may forward a warranty claim.

- Claims are to be forwarded as soon as a defect is discovered and not later than 2 weeks after the defect is discovered.
- Claims are to be addressed to the sales agent of the wheelchair. Please note that sales documentation has to be filled in and signed correctly with serial number and eventual number in order to document time and place of the purchase of the wheelchair.

- The sales agent and Loh Medical are to decide whether a defect is covered by this warranty. The claimer will be informed about the decision as soon as possible.
- If the claim is accepted, the sales representative are to decide if the product will be repaired, replaced or if the customer is entitled to a reduced price.
- If a warranty claim is judged to be invalid - after careful inspection of the defect (defect due to wrong use and/ or lack of required maintenance) you are free to decide if you want to have the product repaired (if possible) at your expense, or if you want to purchase a new product.

Normal wear, incorrect use or incorrect handling is not a reason for claims.

11.4 Customized / individual adaptations

Coyote T+ Customized / individual adaptations are defined as all adjustments that are not included in this manual. Individual adaptations made by LohMedical are labelled with a unique SN number for identification.

Wheelchairs that are especially adjusted / adapted by the customer cannot keep the CE mark given by Loh Medical. If the adjustments are performed by other than Loh Medical approved dealers, the warranty given by Loh Medical will not be valid.

If you have any uncertainty about special fitting and adaptations, please contact Loh Medical.

If you have different needs than what our standard wheelchair program can cover, please contact customer service for eventually special adjustments or solutions.

11.5 Combinations with other products

Combinations of **Coyote T+** and other products not manufactured by Loh Medical: Generally in these cases, the CE mark of all the products involved will not be valid. However, Loh Medical has made combination agreements with some manufacturers about some combinations. By these combinations the CE mark and warranties are valid.

For further information, please contact your dealer or Loh Medical directly.

PRODUCT RESPONSIBILITY

Coyote T+ with different configurations of equipment has been tested and risk evaluated by Loh Medical.

Any alterations or substitutions must not be made to the wheelchair securement points or to structural and frame parts without consulting the wheelchair manufacturer Loh Medical.

-  Substitutions or alterations of components from third part suppliers to **Coyote T+** requires the risk evaluation and acceptance of the product responsibility and safety for use of the wheelchair from the manufacturer that is performing the substitution or alteration.

11.6 Service and repair

Information about service and repair services in you area, please contact your local dealer.

-  **A unique identification number / serial number is to be found on the crossbar on the wheel frame on left side of the chair.**
-  **A spare part catalogue for the wheelchair can be obtained through your local dealer or downloaded at www.lohmedical.com**
-  **A refurbishment manual for the wheelchair can be obtained through your local dealer or downloaded at www.lohmedical.com**
-  **Information on product safety notices and product recalls are available at www.lohmedical.com**
-  **A recycling manual for the wheelchair can be obtained through your local dealer or downloaded at www.lohmedical.com**

12. MEASUREMENTS & WEIGHT

Size - SW*	Size depth**	Back height *** (Extender)	Total width**** (+230)	Weight ex cushion
350 mm	(375) 400 - 500 mm	515 (625) mm	580 mm	34 kg
400 mm	(375) 400 - 500 mm	515 (625) mm	630 mm	35 kg
450 mm	(375) 400 - 500 mm	515 (625) mm	680 mm	36 kg
500 mm	(375) 400 - 500 mm	515 (625) mm	730 mm	37 kg

* SW - Seat Width: Dimension between edge of frame tubes.

The distance between cloths protectors: remove 15 mm.

** Dimension from front edge of the seat plate to the back hinge – without cushion.
By correct placed UNO back cushion ca. 30 mm must be subtracted from this dimension.

*** Dimension from seat plate to top of the back.

**** Minimal total width = SW + 190 mm. Max total width without other adj. = SW+230mm.

With adjusted SW (+50mm), the total width = SW+290.

 **Weight include main wheels, front castors, leg supports and arm supports.**

 **Recommended tyre pressure for pneumatic drive wheel tyres: 40–45 PSI.**

 **Max. user weight for Coyote T+ with seat width up to 500 mm 136 kg.
When used as seat in a car: Max. user weight 136 kg.**

 **When accessories and extra equipment are mounted, the weight of them must be subtracted from the max user weight.**

 **Luggage loaded on to the wheelchair must not exceed 10 kg.
The luggage must not be placed in a manner that reduce the stability of the chair.**

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