

Zippie Zone March 2018

FIR71

	<u> </u>				LINZ	.,					
	Date:							SUBMIT	TING FOR:	□ Quote □	Order
	Dealer Acct #:							OOD		_ Q0010	Ordor
	Dealer:							PO#:			
_	ealer Contact:							го# p To:			
D	ealer Address:		0.7					ntion:			
	Dealer City:	, ,	ST: ZIP:					ress:			
	Dealer Phone:	()	Fax: ()					ress:			
							Ship To	· —		ST:	ZIP:
	Confirm Via:	□ Fax	☐ Email				Ship To Ph	none: ()	Fax: ()	
	Mark For:										
			The HCPCS CODES provided are based	on code verificati	on thro	uah P	DAC or our inte	erpretation	of the code definit	tions.	
		TI	he use of HCPCS CODES does not ensure covera								
~	Option #	HCPCS	Definition	Price		~	Option #	HCPCS		Definition	Price
•	opaon n	7.0.00	Beimaen	1. FRAN	/IE ST	YLE	option ii	7.07.00		Dem more	7 1100
			Frame w/ seating		1 [Frame w/	o seating	
	190FS38	E1235	Zippie Zone w/ seating	\$2,515			190FS40	E1237	Zippie Zone w/	o seating	\$2,365
			Zippie Zone, Flared w/ seating		-				Zinnie Zone E	lared w/o seating	
	190FS39	E1235	(Increase of total frame length appx. 2" - 1" per	\$2,515			190FS41	E1237		I frame length appx. 2" - 1" per	\$2,365
			side flare.)	. ,					side flare.)		. ,
	Note:	: HCPCS Code	es shown as submitted. Some sizes fit other HCPC	;	-	-	Note: I	HCPCS Co	des shown as sub	omitted. Some sizes fit other HCP	'n
		C	codes i.e. 15"-16" widths = K0005						codes i.e. 15"-16	" widths = K0005	
				0.75		-					
				2. 18	ANSI	1					
			Transit Option								
	190TO41	K0108	(Rear frame extends 1" to accom. transit & not lose depth adj. range; N/A w/ 18" depth; N/A w/	\$300							
			seat ht. below 15" w/ inset frame)								
			,		-						
			3. FRAME SETUP/ ANGLE		1 F				4. ZIPPIE	PACKAGE	
	190FA8		80° Frame Angle	NC	J L		190QE01		Z Performance		
_	1001710		or maineraligie		-	_	.000201		2 : 0::0:::::a::00	Single Sided Fork	\$125
										Ultra Lightweight Wheel	\$250
		Fr	ame Inset Per Side (Footrest)		1					Spoke Guards Clear Zippie	\$115
	190FF1		0"	NC	_					Compact Wheel Lock	\$125
			1"		-					Retail Value	\$615
	190FF2		(N/A 8"-9" frame width or with flared frame)	NC						Package Price	\$310
				3. FRAM	E CO	L OB					
				J. FRAIN	IE COI	LON					
			Standard Colors		J				Standard M	latte Colors	
	190C140		Aztec Gold	NC	_		190C116		Matte Black		NC
	190C118		Bead Blast	NC	_		190C151		Matte Black Ch	nerry	NC
	190C2		Black (N/A as base w/ Zebra)	NC	_		190C152		Matte Purple		NC
	190C141		Black Cherry	NC	_		190C150		Matte Electric	Blue	NC
	190C10		Black Opal	NC	_		190C149		Matte Evergree	en	NC
	190C38		Blue Opal	NC							
	190C36		Candy Blue	NC					Custom	n Colors	
	190C14		Candy Purple	NC			190C43		Mossy Oak Ca	mouflage	\$295
	190C35		Candy Red	NC			190C42		Stars and Strip	es	\$295
	190C142		Electric Blue	NC		_	1000109		Zebra (Choose	base color above: White Std.	\$20E
	190C143		Evergreen	NC	_		190C108		N/A black)		\$295
	190C126		Glow	NC	_						
	190C144		Green Apple	NC					Anodized Co	olor Package	
	190C147		Hot Sparkle Pink	NC	٠.		Forks, Backres	t Bracket, I		If Frog Legs are chosen, forks w	ill be black
	190C148		Mauve Pink	NC	_		190EB07		Black		NC
	190C103		Rootbeer	NC	_		190EB08		Blue		\$150
	190C114		Sparkle Silver	NC	_		190EB13		Gold		\$150
	190C145		Sunrise Orange	NC	-		190EB09		Green		\$150
	190C113		Titanium Color Paint	NC	-	\Box	190EB11		Red		\$150
	190C8		Yellow	NC	-		190EB12		Ti Grey		\$150
					-	-					

	4 FRANCIUM TO THE TOTAL THE TOTAL TO THE TOTAL THE TOTAL TO THE TOTAL THE TOTAL TO THE TOTAL TOT												
	4. FRAME WIDTH												
	190W35		8" (N/A w/ 17"-18" depth)	NC		190W4		13"	NC				
	190W36		9" (N/A w/ 17"-18" depth)	NC		190W5		14"	NC				
	190W1		10" (N/A w/ 17"-18" depth)	NC		190W6		15"	NC				
	190W2		11" (N/A w/ 17"-18" depth)	NC		190W7		16"	NC				
	□ 190W3 12" NC												
	5. SEAT DEPTH (FRAME DEPTH)												
		For additional	caster/wheel clearance you can choose a differen	t frame length and sl	ing depth	(i.e. shorter sli	ng vs longe	er frame to allow more wheel/caster combos.)					
			8" (Frame grows 9" -11")			190SD7		14" (Frame grows 15" - 17")	NC				
	190SD1		(Only avail. With 18" rear wheels	NC		190SD8		15" (Frame grows 16" - 18")	NC				
			and 3" casters)			190SD9		16" (Frame grows 17" - 18")	NC				
Ξ.	190SD2		9" (Frame grows 10" - 12")	NC NC		190SD10		17" (Frame grows 18")	NC				
	190SD3 190SD4		10" (Frame grows 11" - 13") 11" (Frame grows 12" - 14")	NC NC		190SD11		18" (No Depth Growth) (N/A with Transit)	NC				
<u>.</u>	190SD4		12" (Frame grows 13" - 15")	NC NC				(IVA WIII Transit)					
	190SD5		13" (Frame grows 14" - 16")	NC									
Ц.	190300		13 (Frame grows 14 - 16)										
	6. SLING DEPTH												
		For additional	caster/wheel clearance you can choose a differen			<u> </u>	ng vs longe	, , , , , , , , , , , , , , , , , , ,					
Π.	190GJ1		8" (Sling Grows 9"-10")	NC		190GJ6		13" (Sling Grows 14"-15")	NC				
	190GJ2		9" (Sling Grows 10"-11")	NC		190GJ7		14" (Sling Grows 15"-16")	NC				
	190GJ3		10" (Sling Grows 11"-12")	NC		190GJ8		15" (Sling Grows 16"-17")	NC				
	190GJ4		11" (Sling Grows 12"-13")	NC		190GJ9		16" (Sling Grows 17"-18")	NC				
	190GJ5		12" (Sling Grows 13"-14")	NC		190GJ10		17" (Sling Grows 18")	NC				
						190GJ11		18" (No Growth)	NC				
	7. SEATING												
			Seating Systems					Width					
			Ultra Seat Sling					8" 9" 10" 11" 12" 13" 14" 15" 16"					
	190OP05		(Tension Adjustable)	NC			8" 9"	8" 9"					
	190OP04		Omit Seat Sling	NC			10"	10"					
	190CU30	K0669	1" Cushion	\$120									
		1.0000	(See matrix for availability)	V.20			12" 12"	11"					
	190CU1	K0669	2" Cushion	\$120			13" 14" 15"	13" O D L L L L L L L L L L L L L L L L L L					
			(12" wide and up only)										
			IAV Cuchions (Ontional)		_		16"	16" 17"					
		To order any	JAY Cushions (Optional) IAY cushions on-chair, please see the JAY price lis	.+			17" 18"	17"					
		•	.,					8" 9" 10" 11" 12" 13" 14" 15" 16"					
	Cusn	JAY Bundle	art number must be filled out with desired width and	и аерт	_			Width					
	190JC20		ime or Part #	Price List				All Cooking and lable at the continue					
_	1300020	Width (in.):	Depth (in.)	T HOC EIST				1" Cushion available at these sizes 1" Cushion NOT available at these sizes					
		Widai (iii.).	Борит (пт.)	0. DAOM	DECT								
				8. BACKI	VE91								
		1	Backrest Type					JAY Backs (Optional)					
	190BT96		Folding Angle Adj. w/ Double Lock Down	NC	To c	rder any JAY		-chair, please see the JAY Price Lst. Require	es Part #				
			,		_			ndled Backrest					
						190JB02	Backres	t Name or Part #	Price List				
Muset	aammilata Caatia	1	Backrest Upholstery				Width (in.) Height (in.):					
iviust	complete Sectio	п а оп пехі ра	1-				Paol	krest Handle Option (Stroller Handle)	1				
	190BU1		Standard Upholstery (N/A with Fold Down Push Handles)	NC	8"-10)" widths have a		iller handle; omit push handle recommended. Rece	eiver places				
								n and inside 4" ext. Double "D" tubing to eliminate r					
	190BU72		Omit Std. Upholstery (N/A with Fold Down Push Handles)	NC		190A192	Т.	Telescoping Push Handles (19" from Seat)	\$215				
	1900072		(Default when choose w/o seating)	NO		130/132			Ψ213				
Must	complete Sectio	n B on next na	1,			190A195		Telescoping Push Handles - Tall (25" from seat to push handle)	\$215				
		<u> </u>						,,					
	190BU70	E2611	Tension Adj. Ballistic Nylon 186BU70 (N/A with 8"-9" widths)	\$315				Backrest Extension					
	190BU73		Omit Tens. Adj. Ballistic Nylon	NC	Stan	dard Extension	2"	Daoniosi Extension					
Δ.	1000010		Crime Forto. Aug. Damotto Hyloti	140	Starr	aa. a Exterioren	_	Rear Backrest Extension					
						190OD20		Brackets (4" overall)	\$95				
								2.2.1 2.2.2					

1a: OMIT PUSH HANDLE OPTIONS W/ STANDARD UPHOLSTERY OR WHEN OMITTII	NG UHOLSTERY
□ 190PH2 Omit Push Handles NC	
Choose 0 or 1" Inset and Back Height	
□ 1900D7 0 inset (straight back tubes) NC □ 1900D5 1"	inset (N/A 8"-9" widths) \$80
The first grant and grant	"- 13" Back Height NC
	1 11"
	" - 16" Back Height NC] 14"
Setting 13" - 17" Back Height NC Setting 16	" - 18" Back Height NC
] 16"
□ Setting □ 15" □ 16" □ 17" □ 18"	
2A: INTEGRAL PUSH HANDLE OPTIONS W/ STANDARD UPH.	
□ 190PH6 Integral Push Handles NC	
Choose 0 or 1" Inset and Back Height	
□ 1900D7 0 inset (straight back tubes) NC □ 1900D5 1"	inset (N/A 8"-9" widths) \$80
11 Setting	"- 14" Back Height NC
J L 9" L 10"	1 12" □ 13" □ 14" " - 16" Back Height NC
Setting] 14"
	" - 18" Back Height NC
15" - 18" Back Height NC] 16"
□ Setting □ 15" □ 16" □ 17" □ 18"	
SECTION B - Tension Adjustable Ballistic Nylon	1
4D. CAN'T DUCKLUANDLE OPTIONS W/ TENSION AD L	
1B: OMIT PUSH HANDLE OPTIONS W/ TENSION ADJ.	
□ 190PH2 Omit Push Handles NC	
☐ 190PH2 Omit Push Handles NC Choose 0 or 1" Inset and Back Height	
□ 190PH2 Omit Push Handles NC Choose 0 or 1" Inset and Back Height □ 1900D7 0 inset (straight back tubes) NC □ 1900D5 1"	inset (N/A 8"-11" widths) \$80
□ 190PH2 Omit Push Handles NC Choose 0 or 1" Inset and Back Height □ 1900D7 0 inset (straight back tubes) NC □ 1900D5 1" □ Setting 11"-14" Back Height NC □ Setting 11	inset (N/A 8"-11" widths) \$80 "- 13" Back Height NC 1 11"
□ 190PH2	"- 13" Back Height NC 1 11"
190PH2	"- 13" Back Height NC 1 11"
□ 190PH2	"- 13" Back Height NC 1 11"
□ 190PH2	"- 13" Back Height NC 1 11"
190PH2	"- 13" Back Height NC 1 11"
190PH2	"- 13" Back Height NC 1 11"
□ 190PH2 Omit Push Handles NC	"- 13" Back Height NC 1 11"
190PH2	"-13" Back Height NC 1 11"
190PH2	"- 13" Back Height NC 1 11"
190PH2	"-13" Back Height NC 1 11"
190PH2	"-13" Back Height NC 11"
190PH2	"-13" Back Height NC 11"
190PH2	"-13" Back Height NC 11"
190PH2	"-13" Back Height NC 11"
190PH2	"-13" Back Height NC 11"
190PH2	"-13" Back Height NC 11"
190PH2	"-13" Back Height NC 11"
190PH2	"-13" Back Height NC 1 11"
190PH2	"- 13" Back Height NC
190PH2	"-13" Back Height NC
190PH2	"- 13" Back Height NC

SECTION A - Standard Upholstery

_													
							9.	ARMREST	S				
		A	rmrest - Pa	dded Swin	ig Away						Armrest	- T-Post Foam (see armrest matrix below)	
				e 1" Seat De						190AA129		T-Post Foam xtra tall	\$210
	190AA23			wing Away			\$210	<u> </u>		190AA130	E0973	T-Post Foam 2x tall	\$210
	[Setting	□ 5"	□ 6"	□ 7"	□ 8"	□ 9"						
	190AA18	E0973	Doddod C	wing Away	Ctd		\$210	L				Height Adjustable (see armrest matrix below) std. armrest receiver only - N/A with 8"-9" depths	
	190AA 16	E0973	□ 9"	wing Away	- Siu.	□ 12"	□ 13"				onies with		
	[☐ Setting	□ 14"	<u> </u>	<u> </u>	□ 1Z	L 13			190AA11	E0973	Single-Post Height Adj. Std. (N/A with 8"-11" depths)	\$285
	Wit	h Transit, Lose 1	" Seat Depth	Growth, belo	ow option add	ds 3" in Leng	gth			190AM23		Desk Classic	NC
	190AA127	7 E0973	Padded S	wing Away	+3" - Low		\$210)	\Box	190AM21		Desk Contour	NC
_	[Setting	□ 5"	□ 6"	□ 7"	□ 8"	□ 9"			190AM24		Full Classic	NC
	190AA128	B E0973	Padded S	wing Away-	+3" - Std.		\$210			190AM22		Full Contour	NC
	[Setting	□ 9" □ 14"	□ 10"	□ 11"	□ 12"	□ 13"			190AA29	E0973	Single-Post Height Adj. Low. (N/A with 8"-9" depths)	\$285
-									ο.	190AM23		Desk Classic	NC
										190AM21		Desk Contour	NC
										190AM24		Full Classic	NC
										190AM22		Full Contour	NC
							est Matrix	for Single I	Post	and T-Post			•
Ar	mrest 7	ransfer Tube			CLAS	SSIC PAD				-	L		
	Size	Length (D)		Height Ran	ge (H)		Pad Lei Desk	ngth (L) Full				Single	Post
	STD	8"		8" to 1	3"		10-1/4"	14"					
L	.OW	5-1/4"		6-1/2" to 1	0-1/2"		10-1/4	14					
											1		
Ar	mrest 7	ransfer Tube			CONT	OUR PAD		tl- /l \					
	Size	Length (D)		Height Ran	ge (H)	-	Pad Lei Desk	Full				- D -	
	STD	8"		7-1/2" to 1:	2-1/2"								
L	.OW	5-1/4"		6" to 1	0"		10-1/2"	14-1/4"	'			٥٠	H
			Arm	rest - T-Po	st Foam							افض	
	Armrost Sizo Height Pango (H) Pad Langth (L)				ao (H)		Pad Lo	agth (L)					

7-1/2"

7-1/2"

Xtra Tall

2X Tall

6" to 7"

7" to 8"

T-Post

10. FOOTRESTS

Footrest											
190F228		Aluminum Tubular	NC								
190F273		Aluminum Tubular w/ Snap-On Cover	\$30								
190F50	K0040	Aluminum Platform (8" width is only available with high mount option)	\$95								
190F276	K0040	Angle Adj. Flip Back (N/A with 6" caster with an 8" width)	\$185								

Adaptive Platform Spacers										
N/A with Aluminum Tubular. Now available with Angle Adj. Flip Back.										
190AD1	Adaptive Platform 2" Spacer	\$95								
190AD2	Adaptive Platform 4" Spacer	\$95								
190AD3	Adaptive Platform 6" Spacer	\$95								

	Footrest Options (see matrix below)									
□ 190Dl08			High Mount Ext Tubes (N/A w/ tubular footrest; N/A w/ 16" wide flared frame) **	\$55						
	190DI03	K0038	Leg Strap	\$35						
	190DI04	K0038	Extra Leg Strap	\$35						
	190DI05		Gel Hanger Cover	\$105						
	190DI07		Neoprene Impact Guards	\$85						
	190DI02		Wear Pads	\$50						

SEAT SLING TO FOOTREST DISTANCE (based on 1" Seat Angle)

*Seat to Footrest distances may be +/- 1" from data shown dep. on frame type.

** High mount can adj. full range to seat sling

Note: Flip back may interfere with sling at highest settings.

Note: Flip back may interfere with sling at highest settings.									
FSH	Caster	*Minimum	*Maximum						
14.5*	3-	10.4"	11.4"						
14.5*	4-	9.4*	10.9"						
15.0"	3"	11.4"	12.4"						
15.0"	4"	10.4"	11.9"						
15.0"	5*	9.4"	11.4"						
15.5"	3"	11.4"	12.4"						
15.5"	4"	10.4"	11.9"						
15.5"	5*	10.4"	12.4"						
15.5"	6"	9.4"	11.9"						
16.0"	3"	12.4"	13.4"						
16.0"	4"	12.4"	13.9"						
16.0*	5*	10.4"	12.4"						
16.0*	6"	10.4"	12.9"						
16.5*	3.	12.4"	13.4"						
16.5*	4"	11.4"	12.9"						
16.5*	5"	10.4"	12.4"						
16.5*	6"	10.4"	12.9"						
17.0*	3•	13.4"	14.4"						
17.0*	4"	12.4"	13.9"						
17.0*	5"	11.4*	13.4"						
17.0*	6"	11.4*	13.9"						
17.5*	3"	13.4"	14.4"						
17.5*	4"	12.4"	13.9"						
17.5*	5*	11.4*	13.4"						
17.5*	6"	11.4"	13.9"						
18.0*	3"	14.4"	15.4"						
18.0*	4"	13.4"	14.9"						
18.0*	5*	12.4"	14.4"						
18.0*	6"	11.4"	13.9"						
18.5*	3"	14.4"	15.4"						
18.5*	4-	14.4"	15.9"						
18.5*	5*	12.4"	14.4"						
18.5*	6"	12.4"	14.9"						
19.0*	3•	14.4"	15.4"						
19.0*	4"	14.4"	15.9"						
19.0*	5*	12.4"	14.4"						
19.0*	6"	12.4"	14.9"						
	•		-						

Not Available with 1" Inset

Front Seat Height Matrix	
--------------------------	--

Front Seat Height Matrix Front Caster									
Rear Wheel Size	□ 3"	□ 4"	□ 5"	□ 6"					
	☐ 14.0"	☐ 14.0"	N/A						
	□ 14.5"	□ 14.5"	N/A	N/A					
	□ 15.0"	□ 15.0"	□ 15.0"						
	☐ 15.5"	□ 15.5"	□ 15.5"	□ 15.5"					
	☐ 16.0"	□ 16.0"	□ 16.0"	☐ 16.0"					
18"	□ 16.5"	□ 16.5"	□ 16.5"	☐ 16.5"					
	□ 17.0"	□ 17.0"	□ 17.0"	□ 17.0"					
	□ 17.5"	□ 17.5"	□ 17.5"	□ 17.5"					
	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.0"					
	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 18.5"					
	☐ 19.0"	☐ 19.0"	□ 19.0"	□ 19.0"					
	☐ 14.0"	☐ 14.0"	N/A						
	☐ 14.5" —	☐ 14.5" —		N/A					
	☐ 15.0"	☐ 15.0"	☐ 15.0"						
	☐ 15.5" —	☐ 15.5"	☐ 15.5"	☐ 15.5"					
	☐ 16.0"	☐ 16.0"	☐ 16.0"	☐ 16.0"					
20"	☐ 16.5"	☐ 16.5"	☐ 16.5"	☐ 16.5"					
	☐ 17.0"	☐ 17.0" —	☐ 17.0"	☐ 17.0" —					
	☐ 17.5"	☐ 17.5" —	☐ 17.5"	☐ 17.5" —					
	☐ 18.0"	☐ 18.0"	□ 18.0"	☐ 18.0"					
	☐ 18.5"	□ 18.5"	□ 18.5"	□ 18.5"					
	□ 19.0"	□ 19.0"	☐ 19.0"	□ 19.0"					
	☐ 15.5"	☐ 15.5"	☐ 15.5"	☐ 15.5"					
	☐ 16.0"	☐ 16.0"	☐ 16.0"	☐ 16.0"					
	□ 16.5"	□ 16.5"	☐ 16.5"	□ 16.5"					
22"	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"					
	☐ 17.5"	☐ 17.5"	☐ 17.5"	☐ 17.5"					
	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.0"					
	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 18.5"					
	☐ 19.0"	☐ 19.0"	☐ 19.0"	☐ 19.0"					
	☐ 16.0"	☐ 16.0"	☐ 16.0"	☐ 16.0"					
	☐ 16.5"	☐ 16.5"	☐ 16.5"	☐ 16.5"					
24"	☐ 17.0"	☐ 17.0"	☐ 17.0"	☐ 17.0"					
24"	☐ 17.5"	☐ 17.5"	☐ 17.5"	☐ 17.5"					
	☐ 18.0"	☐ 18.0"	☐ 18.0"	☐ 18.0"					
	☐ 18.5"	☐ 18.5"	☐ 18.5"	☐ 18.5"					
	□ 19.0"	□ 19.0"	□ 19.0"	□ 19.0"					

FRONT SEAT HEIGHTS ARE SHOWN IN THE BOXES ABOVE.

DEPENDING ON CG SETTING, SEAT DEPTH, AND FRAME SIZE SOME

OF THE ABOVE SIZES MAY NOT BE AVAILABLE

11. CASTERS Caster Wheel 190CW7 5"x1.5" Semi-Pneumatic \$85 190CW1 3" Micro NC 190CW66 5"x1 5" Aluminum Soft Roll \$195 4" Low-Profile Polyurethane 190CW3 NC П 190CW34 6" Polv NC П 190CW31 \$70 190CW35 6"x1.5" Semi-Pneumatic \$85 П 4" Micro Light Up 190CW64 4" x 1.5" Aluminum Soft Roll \$185 190CW5 5" Low-Profile Polyurethane NC Caster Forks and Suspension NC 190CW51 5" Micro Light Up \$85 190CF65 Multiple Position Multi Position Single Sided П 190CF69 \$125 5" Soft Roll (Grey) 190CW6 F2219 \$90 (N/A with Frog Legs) П 190CF53 E1015 Frog Legs 0 - 100 lbs. \$425 П 190CF70 E1015 Frog Legs 101 - 165 lbs. \$425 12. AXLES AND WHEELS Axle Plate Rear Wheels 18" Spoke (N/A with Spinergy, Natural Fit Adjustable 190RW12 \$100 190AL61 NC Handrim or Hub Lock) (COG Adjustable -0.25" through 2.25" in 0.25" increments) 190RW3 20" Mag NC 190RW8 22" Mag NC **Center of Gravity Preset** 190RW1 24" Mag NC Dependent on wheels, casters & depth 190RW103 20" Mag Wheel 5 Spoke NC 190CC21 NC NC -0.25190RW104 22" Mag Wheel 5 Spoke П 190CC26 NC П NC 0.00" (Flush w/ backrest) 190RW102 24" Mag Wheel 5 Spoke 190CC22 0.25 NC 190RW74 NC 20" Lite Spoke Wheel 190CC16 0.50" NC 190RW75 22" Lite Spoke Wheel NC 0.75 NC NC 190CC23 190RW76 24" Lite Spoke Wheel 190CC20 1.00' NC 190RW06 22" Ultra lightweight Spoke Wheel 2.0 \$250 190CC24 1.25" 190RW10 24" Ultra Lightweight Spoke Wheel 2.0 \$250 NC П 190CC15 1.50 NC 190RW28 22" Spinergy SPOX Black Spokes \$820 24" Spinergy SPOX Black Spokes П 190CC25 1.75 NC 190RW95 \$820 NC 190RW32 22" Spinergy LX Black Spokes \$1,025 190CC17 2.00 190CC27 2.25 NC 190RW98 24" Spinergy LX Black Spokes \$1,025 Rear Wheel Tire Camber NC 190RT3 190CI1 00 Full Poly Tire NC 190CI2 30 NC E2211 190RT6 High Pressure Clincher \$95 F2212 6° (N/A with Std Mag, 5 spk MAG Wheels and 190CI4 NC Hub locks or Side guard Hub locks) П 190RT38 E2218 SHOX Performance Poly \$95 F2211 190RT27 Knobby V-TRAK by Kenda \$95 Rear Seat Height E2212 Can not exceed 2.5" drop 190RT1 NC Pneumatic 190RH11 13" NC 190RH12 13.5 NC E2211 190RT2 Pneumatic Airless Insert \$80 F2213 NC 190RH13 14" 190RH14 14.5 NC П E2211 190RT21 Schwalbe Marathon Plus Evolution \$135 E2212 П 190RH15 15" NC 190RH16 NC 15.5 E2211 190RT28 Schwalbe Speed \$125 E2212 190RH17 16" NC 190RH18 16.5 NC 190RT4 Low Poly NC 190RH19 17" NC Whee Ultra 190RH26 NC 17.5 5 Spoke Lite Mag Light-Spinergy Spoke 190RH20 NC Mag Spoke 18" w eight Tire 190RH27 18.5 NC 20",22",24" 20",22",24 20",22 N/A N/A Full Poly Tire 19" 190RH21 NC High Pressure Clincher N/A 24" 22",24 22",24 20",22",24 24" 24" 24" SHOX Performance Poly 24' 24" Knobby V-TRAK N/A N/A **Rear Seat Height Minimum** N/A 20",22",24 20",22",24 20",22",24 Pneumatic 22",24 18' Rear Wheel 14.5" 15.0" 15.5" 16.0" 16.5" 17.0" 18.5" 13.0" 13.5 14.0 17.5 18.0" 19.0

18

20'

22

20",22",24

N/A

N/A

24"

Pneumatic Airless Insert

Schwalbe Marathon Plus

Schwalbe Speed Run

Low Poly

20".22".24"

N/A

N/A

20",22",24"

22" = 489 mm 24" = 540 mm

20".22".24

24"

24"

20",22",24"

22",24

24"

N/A

22",24

22",24"

22".24"

N/A

18'

N/A

^{*} With 22" & 24" rear wheels and KIK or HP Clinchers, you may get an additional 1/2" lower seat height in some scenarios.

12. AXLES AND WHEELS (Cont.)

22" Rear Wheels w/ 5" Casters & 3° or 6° of Camber 0" Frame Inset

Not available in 12" or less seat depths

13" Seat Depth

- Must have front seat to floor height at least 17" in the front and 15-1/2" in the rear.
- You can increase STFs but must keep 1-1/2" dump.
- If you change camber to 0° you can get away with 1" dump
- You must choose COG setting of -1/4".

14" Seat Depth

- Able to get away with 1" seat dump with a minimum of 15-1/2" in the Rear and 16-1/2" in Front.
- Can have COG set as far forward as +1/4".

15" Seat Depth

- Able to get away with 1" seat dump with min of 15" Rear STF and 16" Front STF
- COG can be adjusted as far forward as 1".

22" Rear Wheels w/ 4" Casters, 3° or 6° of Camber, and 0" Frame Inset

Not available in 8" or 9" seat depths.

10" Seat Depth

- Must have at least a 15" Rear STF and 16" Front STF.
- COG setting must be -1/4".
- If you increase Rear STF to 15-1/2" and 16-1/2" Front STF you can move COG to 0".
- If you increase the dump to 1 ½" (15 R-STF and 16 1/2" F-STF) or change Camber to 0 degress you can increase COG to 1/4".
- 11" Seat Depth
- Same as 10" Seat Depth, just allows you to have COG setting a 1/4"

12" Seat Depth

- Must have at least 15" Rear STF and 15-1/2" Front STF.
- COG setting can be as far forward at 1/2".
- If you increase the dump to 1" instead of 1/2" you can increase COG as far forward as 1".

24" Rear Wheels w/ 4" Casters, 3° or 6° Camber, and 0" Frame Inset

Not available in 10" seat depths or less.

11" Seat Depth

- Must have at least a 16" Rear STF and 17" Front STF.
- COG setting can be 0" or -1/4".

12" Seat Depth

- Must have at least 16" Rear STF and 17" Front STF.
- COG setting can be as far forward as 1". (Camber must be 3 degrees)
- If camber is 6° than COG can only be as far forward as 3/4"

	Handrim									
	☐ 190HR1 Aluminum Anodized									
	190HR4	Plastic Coated (N/A on Ultra)	\$95							
	190HR22	Natural Fit -Standard Grip (N/A on Ultra Lite, 20" Lite Spoke or 20" and 22" 5 spk MAG)	\$395							
_	190HR23	Natural Fit - Super Grip (N/A on Ultra Lite, 20" Lite Spoke or 20" and 22" 5 spk MAG)	\$425							
	190HR30	Natural Fit LT - Standard Grip (N/A on Ultra Lite, 20" Lite Spoke or 20" and 22" 5 spk MAG)	\$415							
_	190HR31	Natural Fit LT - Super Grip (N/A on Ultra Lite, 20" Lite Spoke or 20" and 22" 5 spk MAG)	\$475							
	190HR24	Spinergy Flex Rim (Only on 24" Spinergy LX)	\$650							
	190HR7	Omit Handrim (Only on Mag, 24 Spinergy)	NC							
		Axle								
	190AX1	Quick Release Steel	NC							

Hub Lock									
Only available with Lite Spokes; N/A with 6° camber									
190LK03	Seat Frame Mount	\$100							
	(Side Guard - See Below)								

	Side Guard (N/A with Ht. Adj. Armrest or T-Post)						
	190SG33	Plastic Standard (N/A 8"-9" depth; N/A w/ low padded swing aways.)	\$140				
	190SG29	Plastic Low (N/A 8"-9" depth; N/A w/ Padded Swing Aways low setting 5&6)	\$140				
	190SG7	Aluminum Folding	\$185				
	190SG37	Aluminum w/ Fender & Hub Lock (Pull to lock - single hand R-side mount) (N/A with seat frame mount hub lock)	\$275				
_	190SG38	Aluminum w/ Fender & Hub Lock (Pull to lock - single hand L-side mount) (N/A with seat frame mount hub lock)	\$275				

Aluminum Side Guard Color Option						
Only available with Aluminum Side Guard						
□ 190GE01 Black						
190GE02		Matching Frame - Std. Colors	NC			
190GE03		Matching Frame - Custom Color	\$50			

Anti-Tip					
190AT1	E0971	Rear Anti-Tips	\$140		

Dynamic Center Rear Caster						
	1908A20	E0971	Dynamic Center Rear Caster	\$295		

Wheel Lock					
	190WL1	High Mount Push	NC		
	190WL2	High Mount Pull	NC		
	190WL30	Compact	\$125		
	190WL10	Omit Wheel Lock	NC		

Extension Handle					
190EH1	E0961	6" Pair	\$75		

13. ACCESSORIES

Back Pack (Kids Backpack included)					
	190A1	Utility Backpack - Black	\$65		
	190A2	Additional Kids Black Backpack	\$60		
	190A8	Black Seat Pouch	\$50		

Miscellaneous						
190A203	E2207	Crutch Holder	\$195			
190A83		Tool Kit	\$45			

Positioning Belt (w/ clamps)					
	190A112	E0978	1" Dual Pull Side Squeeze Padded	\$85	
	190A117	E0978	Belt 2" Aircraft Buckle	\$110	
	190A118	E0978	Belt 2" Aircraft Buckle Padded	\$120	
	190A141	E0978	Belt 1.5" Auto Buckle w/ clamps	\$70	
	190A814	E0978	1" Comfort Fit APB5116	\$120	
	190A120	K0108	Positioning Belt Side Clamp	\$50	

Spoke Guard						
Not	Not available with 5 spoke MAG wheels					
	190A23	K0065	Black	\$115		
	190A303	K0065	Spoke Guards Clear Zippie (N/A with Spinergy)	\$115		

Upper Extremity Support Tray						
	190A29	E0950	X-Sm (Suggested 10"-12" Width)	\$255		
	190A30	E0950	Sm (Suggested 13"-14" Width)	\$255		
	190A31	E0950	Med (Suggested 15"-17" Width)	\$255		
	190A32	E0950	Lrg (Suggested 18"-20" Width)	\$255		

Additional Options Through Built-4-Me

To inquire about any unique customization for Quickie products, contact: Phone: (800) 333 - 4000, Fax: (800) 333 9011



Commonly Requested Options Include:

- ✓ Aluminum Seat Pan
- ✓ Swing Away Knee Aductor
- ✓ Height Adjustable Stroller Handle
- ✓ Removable Stroller Handle
- ✓ Custom Paint and Anodizing



Sunrise Medical (US) LLC 2842 Business Park Ave. · Fresno, CA 93727 · USA