

# GRAB BARS & ANCHORS



It is a collection of suction cup grab bars with different types of handles and arm supports to assist in the stabilization of the upper extremities of people with disabilities. They increase upper body control to make tasks and activities of daily living easier.

The padded arm Anchor stabilizes the forearm. It can be adjusted lengthwise and the handgrip rotates to any angle. **Image 1.**

The elbow Anchor provides stability to the forearm while leaving the hands free for functional activities. **Image 2.**

This ergonomic hand Anchor can be attached to any smooth surface to achieve stability and positioning. **Image 3.**

This horizontal grab bar can be attached to any smooth surface to assist in functional activities. **Image 4.**



Image 1



Image 2




Image 3




Image 4




# Adapted equipment. **Synonym of service.**

Equipment manufactured and adapted for children and adults with disabilities.

## SPECIFICATIONS

|                    |                         |
|--------------------|-------------------------|
| Hand Anchor        | 1 - K821                |
| Wrist Anchor       | 2 - K822                |
| Horizontal Anchor  | 3 - K823                |
| Right Arm Anchor   | 4 - K811 (S) / K813 (L) |
| Left Arm Anchor    | 5 - K812 (S) / K814 (L) |
| Right Elbow Anchor | 6 - K815                |
| Left Elbow Anchor  | 7 - K816                |
| Grab Bar           | 8 - K824                |

Required user  
measurement



## BENEFITS

They can be used on any table or tray and give the client an anchor to control movements

## Discover our products **ONE**



### DESIGNED FOR LATIN AMERICA

Products designed for the needs of Latin America and the Caribbean, backed by our professional support.



### FUNCTIONAL VERSALITY

Our product line includes wheelchairs, mobility scooters, and lift chairs. Ask your local representative.

**For more information on obtaining this product**  
**Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)

