

# Club Sport

**RGK**



The sports wheelchair that allows you to get active without a major commitment. Only the best materials are used in your Club Sport. Chrome Moly is famous for its strength, durability and is surprisingly not that heavy when worked with properly.

Configurable sizes. **Image 1.**

Suitable for multisport use - e.g. Basketball, Tennis, Rugby. **Image 2.**

Adjustable height backrest and centre of gravity. **Image 3.**

Quick release single anti-tip. **Image 4.**



Image 1



Image 2



Image 3



Image 4



# Keeping you in mind. More options.

Quality products with the required accessories at a great price.


## SPECIFICATIONS

Seat width	13" to 17"
Seat depth	12" to 17"
Frame height	19" to 21" (front)
Frame height	16" to 21" (rear)
Product weight	From 10 kg. / 22 lbs.
Weight capacity	125 kg. / 275 lbs.
Gravity center	Adjustable
Camber	18°
Casters	3"
Accessories	Optional
Available colours	23 available

Regulatory measures



Rear wheels  
Quick release



## BENEFITS

For active people mainly motor injuries.

## Discover our products



### DESIGNED FOR LATIN AMERICA

Products designed for the needs of Latin America and the Caribbean, backed by our professional support.



### FUNCTIONAL VERSALITY

Our product line includes wheelchairs, mobility scooters, and lift chairs. Ask your local representative.

**For more information on obtaining this product**  
**Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

