

Comfortable Backrest



The Comfortable stabilizing backrest provides comfortable support and relief at the same time. The modifiable shape offers personalized adjustment to the user's needs, correcting postural defects and thus improving the control of the user's body.

The materials are certified non-flammable, antibacterial and latex-free.
Image 1.

Three sizes of width and height.
Image 2.

Includes manual air pump, and other optional accessories.
Image 3.

The anatomical shape of the backrest and light weight allow for easy installation in a wheelchair, rehabilitation seat, chair or car.
Image 4.




Image 1




Image 2



Image 3



Image 4



For you. The best technology.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Width	Small 52 cm
	Medium 56 cm
	Large 70 cm

Height	Small 40 cm
	Medium 50 cm
	Large 50 cm

User Weight capacity	100 kg. / 200 lbs.
----------------------	--------------------

Product weight	1.5 kg. / 3.3 lbs.
----------------	--------------------

Accessories	Optional
-------------	----------

Flexible tube for greater stability and rigidity



BENEFITS

Designed for users who require moderate positioning and stability.

Discover our products



DESIGNED FOR LATIN AMERICA

Products designed for the needs of Latin America and the Caribbean, backed by our professional support.



FUNCTIONAL VERSALITY

Our product line includes wheelchairs, mobility scooters, and lift chairs. Ask your local representative.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

