

# Elite Sports Wheelchair

**RGK**



Designed entirely around your measurements and your style of play, the Elite gives you the ultimate responsive connection, ergonomic positioning, and speed that your best performance requires.

The Tri-Frame is a new design that better distributes the energy expended during play. **Image 1.**

Carbon Fiber Sideguards. **Image 2.**

Carbon Fiber Axle. **Image 3.**

Elite CX: Superior Ergonomics & Carbon Fiber Seating System. **Image 4.**




Image 1



Image 2



Image 3



Image 4




# Keeping you in mind. More options.

Quality products with the required accessories at a great price.

## SPECIFICATIONS

Seat width	10" - 20"
Seat depth	10" - 20"
Frame height	14" - 24.8" (front)
Frame height	14" - 24.8" (rear)
Chair weight	8.1 kg.
Weight capacity	125 kg. / 275 lbs.
Center of gravity	Adjustable
Camber	16°, 18°, 20° & 22°
Casters	3"
Rear wheels	24", 25", 26" & 700
Accessories	Optional
Available colors	30 available

Regulatory measures



Rear wheels  
Quick release

## BENEFITS

For active people mainly motor injuries.

## Discover our products



### DESIGNED FOR LATIN AMERICA

Products designed for the needs of Latin America and the Caribbean, backed by our professional support.



### FUNCTIONAL VERSALITY

Our product line includes wheelchairs, mobility scooters, and lift chairs. Ask your local representative.

**For more information on obtaining this product**  
**Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

