

# MaxLite

**RGK**



Only the best materials are used in your MaxLite. Titanium is famous for having all of the positives of aluminium with the additional bonus of a natural shock absorption by the material itself. The utmost best performance of your chair is ensured by only using elements produced by market leaders, alongside a staggering 19 quality checks throughout the build, from measure to handover.

The wheelchair frame is made from commercially pure Titanium. **Image 1.**

With over 30 measurements, can be perfectly tailored to suit the user's disability and clinical requirements. **Image 2.**

Choose an accent colour to suit your style. **Image 3.**

Airtech upholstery system. **Image 4.**



Image 1



Image 2



Image 3



Image 4

# Keeping you in mind. More options.

Quality products with the required accessories at a great price.

## SPECIFICATIONS

Frame Material	Titanium
Frame style	Rigid
Seat Width	10" - 17" / 25 cm - 44 cm
Seat Depth	10" - 17" / 25 cm - 44 cm
User Weight Capacity	100 kg. / 220 lbs.
Product Weight	From 7.5 kg. / 16.5 lbs.
Accessories	Optional
Available colours	25 options



Rear wheels  
Quick release

## BENEFITS

For active people mainly motor injuries.

## Discover our products



### DESIGNED FOR LATIN AMERICA

Products designed for the needs of Latin America and the Caribbean, backed by our professional support.



### FUNCTIONAL VERSALITY

Our product line includes wheelchairs, mobility scooters, and lift chairs. Ask your local representative.

**For more information on obtaining this product**  
**Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

