

# PAL SEATING SYSTEM

**LECKEY**



Reinforced chair that is easy to use at school, home or rehabilitation centers. For mild or moderate postural needs, it improves sitting stability, which allows attention periods to be prolonged and facilitate manipulation activities and fine motor activity.

Two-point lap belt and side pads for lateral stability. **Image 1.**

Height adjustable armrests, legs and footplate. **Image 2.**

Base and footrest options available, including castors, rockers and stabilisers. **Image 3.**

Adjustments Made Easier Than A-B-C. **Image 4.**



Image 1



Image 2



Image 3



Image 4



# Adapted equipment. **Synonym of service.**

Equipment manufactured and adapted for children and adults with disabilities.


## SPECIFICATIONS

|                  |                                    |
|------------------|------------------------------------|
| Sizes            | 4                                  |
| Aproximated age  | 1 to 12 years                      |
| Max user weight  | 30 to 50 kg. /66 up to 110 lbs.    |
| User height      | 75 to 160 cm / 29.5" up to 63"     |
| Seat width       | 15 up to 38 cm / 6" up to 15"      |
| deep seat        | 15 up to 40 cm / 15" up to 15.7"   |
| Backrest height  | 26.5 up to 45 cm / 10" up to 17.7" |
| Backrest recline | From -95° to +105°                 |

Colors



Required user  
measurement  
(total height and  
knee to heel)



## BENEFITS

Pediatric postural seating system, with mild to moderate postural needs.

## Discover our products **ONE**



### DESIGNED FOR LATIN AMERICA

Products designed for the needs of Latin America and the Caribbean, backed by our professional support.



### FUNCTIONAL VERSALITY

Our product line includes wheelchairs, mobility scooters, and lift chairs. Ask your local representative.

**For more information on obtaining this product**  
**Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)

