

Advance Lift



 Joerns

The Hoyer® Advance® is compact, yet able to perform an outstanding range of transfers. Whether from the floor or onto a healthcare bed, the Advance copes with ease. The Advance's unique tapered leg design allows the lift to get even closer to the widest of obstacles. This is especially useful when encountering large chairs, bulky commodes, and wheelchairs. The tapered design also ensures that the resident feels safe and secure during the transfer.

Compact triangular folding design for easy transport and storage. **Image 1.**

Ergonomic design considerations and interaction points. **Image 2.**

Integrated safety system to help protect patients. **Image 3.**

Angled leg design ensures excellent access around furniture. **Image 4.**



Image 1



Image 2



Image 3




Image 4

Adapted equipment. **Synonym of service.**

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS

Closed base	67 cm / 26.3"
open base	108 cm / 42.5"
Maximum height	186 cm / 73"
Minimun height	136 cm / 53.5"
base height	3 cm / 1.1"
Wheel size	3"
Weight capacity	155 Kg. / 341 lbs.
Product weight	31.7 Kg. / 69.9 lbs.
System	Electric (optional hydraulic)
Accessories	Transfer meshes



BENEFITS

Essential in transfers for patients and caregivers in a comfortable and safe way.

Discover our products **ONE**



DESIGNED FOR LATIN AMERICA

Products designed for the needs of Latin America and the Caribbean, backed by our professional support.



FUNCTIONAL VERSALITY

Our product line includes wheelchairs, mobility scooters, and lift chairs. Ask your local representative.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com

