

BEACH ACCESS CHAIR

This sturdy, light-weight, all-terrain, specialized beach wheelchair allows easier access in sand, snow, and other soft soils. Four balloon flotation tires give the chair extra stability that resists tipping and corners better than three wheel designs. Comes standard with single motion park brakes, safety belt, umbrella, sliding footrest, storage bag and comfortable cushion seat. Quickly disassembles into three pieces for easy transport.



Adapted equipment. **Synonym of service.**

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS

Material	Non-corrosive PVC
Wheels	4 wheels 19"
Seat width	20"
Seat depth	20"
Weight capacity	113 kg. / 250 lbs.
Armrests	Folding
Wheel material	POLYVINYL CHLORIDE
Hardware	Stainless steel
Accessory	Seat belt and umbrella

BENEFITS

Excellent help for patients with different transfer needs and limited mobility.

Discover our products



DESIGNED FOR LATIN AMERICA

Products designed for the needs of Latin America and the Caribbean, backed by our professional support.



FUNCTIONAL VERSALITY

Our product line includes wheelchairs, mobility scooters, and lift chairs. Ask your local representative.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

