

Deluxe Lift

 Joerns



Designed to meet the needs of the caregiver as well as the patient, the HPL402 features ergonomic handle grips to reduce back strain and a mechanical down-release lever for non-powered lowering. The single-speed actuator is reverse-mounted to give the patient increased leg clearance, and the unit also features leg guards and caster bumpers to prevent damage to walls and furniture.

Adjustable U-base for easy sitting or toileting and to expand the base with the pull of a lever. **Image 1.**

Battery includes emergency stop button. **Image 2.**

Allows 4-point, 6-point, or chain style slings. **Image 3.**

Battery Type 2-12V, Rechargeable. Charger Input 120V AC. **Image 4.**



Image 1



Image 2



Image 3



Image 4

Adapted equipment. **Synonym of service.**

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS

Closed base	61 cm / 24"
open base	100 cm / 39.5"
Maximum height	198 cm / 78"
Minimun height	77.5 cm / 30.5"
base height	5.5 cm / 2.1"
Wheel size	3"
Weight capacity	181 kg. / 400 lbs.
Product weight	42.2 Kg. / 93 lbs.
System	Electric
Accessories	Transfer meshes



BENEFITS

Essential in transfers for patients and caregivers in a comfortable and safe way.

Discover our products **ONE**



DESIGNED FOR LATIN AMERICA

Products designed for the needs of Latin America and the Caribbean, backed by our professional support.



FUNCTIONAL VERSALITY

Our product line includes wheelchairs, mobility scooters, and lift chairs. Ask your local representative.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

