

HIPPOCAMPE ALL-TERRAIN CHAIR



Since 2002, the Hippocampe wheelchair offers individuals with disabilities accessibility to difficult terrains, impossible in a standard wheelchair: beach, sand, swimming area, hiking trails, snow, mountain.



Adapted equipment. **Synonym of service.**

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS

Frame	Stainless steel and anodized aluminum
Wheels	Optional for pool or beach
Seat width	39.4cm/15.5"
Seat depth	36.8cm/14.5"
Backrest height	36.8cm/14.5"
Total depth	72.4cm/28.5"
Weight capacity	130 kg / 286 lbs.
Accessories	Optional
Frame colors	Blue and pink



BENEFITS

Excellent help for patients with different transfer needs and limited mobility.

Discover our products **ONE**



DESIGNED FOR LATIN AMERICA

Products designed for the needs of Latin America and the Caribbean, backed by our professional support.



FUNCTIONAL VERSALITY

Our product line includes wheelchairs, mobility scooters, and lift chairs. Ask your local representative.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com

