

Hydraulic Patient Lift



The Hydraulic Patient Lift from ProBasics is economically priced and has a 450 lb. weight capacity. The dual purpose padded swivel bar can use 4 or 6 point slings. The rotating handle has a manual hydraulic pump and the frame is a chip resistant brown vein finish.

Heavy-duty round tube steel construction with a scratch-resistant brown finish. **Image 1.**

Made with water-resistant materials and easy to clean. **Image 2.**

Heavy duty 5" solid wheels (2) with reinforced locking. **Image 3.**

Various sizes of nylon four-point strap washable and can be used in wet or dry environments. **Image 4.**




Image 1



Image 2



Image 3



Image 4



Adapted equipment. **Synonym of service.**

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS

ITEM	PBPLSM001
Closed base	56 cm / 22"
Open base	108 cm / 42.25"
Maximum height	162 cm / 64"
Minimun height	92 cm / 22"
Crane height	198 cm / 78"
Wheel size	5"
Weight capacity	204 Kg. / 450 lbs.
Product weight	31 Kg. / 68 lbs.
System	Hydraulic manual
Accessories	Transfer meshes

BENEFITS

Essential in transfers for patients and caregivers in a comfortable and safe way.

Discover our products



DESIGNED FOR LATIN AMERICA

Products designed for the needs of Latin America and the Caribbean, backed by our professional support.



FUNCTIONAL VERSALITY

Our product line includes wheelchairs, mobility scooters, and lift chairs. Ask your local representative.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

