

R20

Empulse



Whether uphill, downhill or on long trips. Foldable and lightweight electric assist device is easily transported, either connected or disconnected to the wheelchair with its practical backpack.

Smooth transport and folding ability.
Image 1.

Simple one-button on/off control for safe handling and easy operation.
Image 2.

Power-assist device enables you to easily ascend slopes of up to 10° with strong 250 W motor. **Image 3.**

Easy assembly. **Image 4.**



Image 1



Image 2




Image 3



Image 4



For you. The best technology.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

ITEM	EISP9
Seat width	13" to 23"
Seat to floor height	15.5" to 20.5"
Distance range	14 km. / 9 miles
Max speed	3.1 mph / 5 km/h
Battery	5.8-6 Ah / 36 V - Lithium Ion
Weight capacity	170 kg. / 375 lbs.
Product weight	5 kg. / 9.9 lbs.



BENEFITS

It is compatible with most folding and tilting manual wheelchairs.

Discover our products



DESIGNED FOR LATIN AMERICA

Products designed for the needs of Latin America and the Caribbean, backed by our professional support.



FUNCTIONAL VERSALITY

Our product line includes wheelchairs, mobility scooters, and lift chairs. Ask your local representative.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com

