

Revolution XL



The new and improved Revolution XL™ Pool Lift is now more sleek and versatile than ever before! Featuring the same robust 500 lb weight capacity as its predecessor, now with increased water draft and anchor setback, this lift provides access to a wider variety of users and provides additional clearance for short pool walls and curbs. The Revolution XL™ is designed to be self-operated by the individual user. With an easy to use 4 button handset, the user is able to control all movement and enjoy a fully independent aquatic experience.

Standard anchor. **Image 1.**

Semi-recessed anchoring. **Image 2.**

Pavement anchor. **Image 3.**

Wood Deck Anchor. **Image 4.**



Image 1



Image 2



Image 3




Image 4



Adapted equipment. **Synonym of service.**

The leaders in recreation, fitness and accessibility, safe and guaranteed.

SPECIFICATIONS

Frame Finished stainless steel
in powder paint

Seat width 55.8 cm / 22"

Seat depth 43 cm / 17"

Backrest height 45.7 cm / 18"

Arm support Folding

foot supports Adjustable and belt

Weight capacity 226 Kg. / 500 lbs.

Rotation 360°

System Electric

Accessories Optional

Colors



Solar charger
helps extend
battery life



Outdoor cover

Match your branding
or decoration with
custom color



BENEFITS

Manufactured for pools and spas, they are ideal for both commercial and residential use.

Discover our products **ONE**



DESIGNED FOR LATIN AMERICA

Products designed for the needs
of Latin America and the Caribbean,
backed by our professional support.



FUNCTIONAL VERSALITY

Our product line includes wheelchairs,
mobility scooters, and lift chairs.
Ask your local representative.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com

