


SUITCASE Ramp



This portable ramp, with its single-fold design, sets up quickly and is ideal for scooters and wheelchairs. Made of lightweight aluminum with an easy-to-carry design, the SUITCASE® Singlefold Ramp is fitted with ergonomically-designed handles and a full-length hinge for additional strength and stability. Designed for high-traffic use, this durable ramp is available in 2' to 8' lengths and features an extruded slip-resistant tread for superior traction.



Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

Material	Aluminum
Broad	75cm/29.5"
Deep	61 - 243 cm / 2' - 8' Ft.
Product weight	3.6 - 13 Kg. / 12 - 30 lbs.
Weight capacity	362 kg / 800 lbs.
Surface	Anti-slip
Creased	In two
Extension	Optional



BENEFITS

Suitable for home, vehicle or travel and offer temporary or long-term solutions

Discover our products



DESIGNED FOR LATIN AMERICA

Products designed for the needs of Latin America and the Caribbean, backed by our professional support.



FUNCTIONAL VERSALITY

Our product line includes wheelchairs, mobility scooters, and lift chairs. Ask your local representative.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

