

Magic 360

magic
Mobility



The Magic Mobility Magic 360 is one of the few compact power chairs that's rugged enough for off-roading, maneuverable enough for indoors, and versatile enough for everything in between.

With a 12" seat lift, the C-Me electric lift function. (5) Travel at 5 km/h with 12" seat elevation. **Image 1.**

Narrow 2" city tires and wide, knobby 4" all-terrain tires provide exceptional grip on all types of different terrain. **Image 2.**

The seat can grow in width from 13" to 23" and grow in depth from 13" to 24", making it perfect for children and adults. **Image 3.**

Offers a wide range of options and accessories. **Image 4.**




Image 1



Image 2



Image 3



Image 4



For you. The best technology.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Seat width	13" up to 23"
Seat depth	13" up to 24"
Product weight	From 145 kg.
Weight capacity	158 kg.
Backrest height	16.5" up to 24"
Maximum rank	33 km.
Maximum speed	8 km/hr.
Wheels	8" and 14"
Turning radius	58 cm (without footrest)
Total width	24" and 26"



- Seat lift and tilt 12" / 45°
- 95-180° Legrest elevation
- Tilt 50°
- 95-160° Recline

BENEFITS

For users with partial or total reduced mobility of lower limbs.

Discover our products



DESIGNED FOR LATIN AMERICA

Our products are designed to meet the needs of Latin America and the Caribbean, offering quality and good prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and motorized wheelchairs, to mobility scooters and standing chairs

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com

