

QUICKIE 2

Q2 Lite / HP

QUICKIE®



With a modular design that can change with you, it offers the widest range of possible configurations in the industry over 100,000 different complex rehabilitation options and counting. Designed to meet unique needs without compromising fit or performance.

Folding frame for easy transport and storage if the user has little space. **Image 1.**

For better positioning and functionality of the chair, it can be configured with JAY backrests and cushions. **Image 2.**

Totally safe quick extraction system on rear wheels for better transportation and assembly. **Image 3.**

Stronger, stiffer crossmember with lighter weight 7000 oval aluminum. **Image 4.**



Image 1



Image 2



Image 3



Image 4




For you. The best technology.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Seat width	12" a 24"
Seat depth	12" a 22"
Product weight	From 9 kg. the frame
Weight capacity	120-158 kg. / 265-350 lbs.
Gravity center	Adjustable 1" - 3.5"
Frame angle	70°, 80°
Foot support	60°, 70°, 80°, Elevating
Accessories	Optional
Settings	Optional

Available colours:



Multiple weight capacity options



Lighter option
QUICKIE 2 Lite

Q-Fit Technology
with perfect fit



Unique
Caster Housing, Fork,
and Wheel Sizes

BENEFITS

For active people mainly spinal and motor injuries.

Discover our products **ONE**



DESIGNED FOR LATIN AMERICA

Our products are designed to meet the needs of Latin America and the Caribbean, offering quality and good prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and motorized wheelchairs, to mobility scooters and standing chairs

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com

