



Mobility & Independence

# QXi



Simplified by Quickie. The Quickie® QXi™ is a durable, low-maintenance frame that is perfect for the person looking for simplicity, longevity, and Quickie's most popular options at the right price.

Durable Frame Design. **Image 1.**

Adjustments, including center of gravity, seat height, camber, and back height. **Image 2.**

Super low seat heights, 13" seat height can be achieved with 5" casters. **Image 3.**

New EXO Upholstery. **Image 4.**



Image 1



Image 2



Image 3




Image 4



# For you. The best technology.

For those who require cutting-edge technology to keep up with the demands of their active life.

## SPECIFICATIONS

Seat width	14" to 20"
Seat depth	14" to 20"
Product weight	From 12 kg. / 26.5 lbs.
User Weight capacity	136 kg. / 300 lbs.
Camber	0° and 3°
Gravity center	Adjustable 1" - 3.5"
Frame angle	60°, 70° and 80°
Foot support	Optional / Elevatable
Accessories	Optional
Settings	Optional

Available colours



## BENEFITS

For active people mainly spinal and motor injuries.

## Discover our products



### DESIGNED FOR LATIN AMERICA

Our products are designed to meet the needs of Latin America and the Caribbean, offering quality and good prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and motorized wheelchairs, to mobility scooters and standing chairs

**For more information on obtaining this product**  
**Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)

