

SR45

QUICKIE®



The simple solution for rotation-in-space technology. The QUICKIE® SR45™ tilt-in-space wheelchair features 45° of Intelligent Rotation in Space Technology and a simplified mix of options to provide proper positioning and comfort at an outstanding value.

45° of intelligent rotation in space technology. **Image 1.**

Durable frame. **Image 2.**

Power tilt, attendant control, toggle switch. **Image 3.**

Quickie's top options. **Image 4.**



Image 1



Image 2



Image 3




Image 4




For you. The best technology.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Frame	Steel
Seat Width	14" to 20"
Seat Depth	14" to 22"
Backrest height	15" to 24"
User Weight capacity	120 kg. / 265 lbs.
Product weight	18 kg. / 39 lbs.
Tilt degrees	Up to 45°
Tilt system	Tilt in space
Accessories	Optional
Settings	In measurements

Available colours:



BENEFITS

A special needs chair that provides tilting and positioning support.

Pedal system
to activate
tilt and brake



Intuitive
foot support designs



Discover our products **ONE**



DESIGNED FOR LATIN AMERICA

Our products are designed to meet the needs of Latin America and the Caribbean, offering quality and good prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and motorized wheelchairs, to mobility scooters and standing chairs

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com

