

Basix2

The lightweight Basix2 wheelchair represents improved adjustability and rigidity, combined with a modern, clean look. Improvements include greater seat depth adjustment (41cm to 46cm), a tension adjustable backrest which is extendable from 41cm to 46cm and a more robust, rigid caster attachment (0, 3 and 6 degrees). Highly efficient and easy to use.



SPECIFICATIONS

- Caster fork adjustment
- Axle plate adjustment
- Elevating footrest
- Rounded front armrest

Color	Silver & White
Material	Aluminium
Max. User Weight	125 kg. / 275.5 lbs.
Seat Angle	0° - 6°
Seat Depths	16", 17", 18"
Seat Height	16.7" - 20" / 42.5 cm - 51 cm
Seat Width	15" - 20.5" / 38 cm - 52 cm
Product Weight	From 14.9 kg. / 32.8 lbs

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com

