

# Economic Folding Chair

 **coyote**



 **ONE**

Coyote is a lightweight, folding wheelchair at a great price. It offers quick-release rear wheels, swing-away armrests, adjustable footplates and an aluminum frame.

Folding aluminum frame for easy transfer and storage if space is limited. **Image 1.**

The armrests are height adjustable to accommodate all users. **Image 2.**

The axle can be adjusted in two dimensions: height and depth (center of gravity). **Image 3.**

The footrest hangers are adjustable in height and can be easily removed. **Image 4.**



Image 1



Image 2



Image 3



Image 4



# Keeping you in mind. More options.

Quality products with the required accessories at a great price.

## SPECIFICATIONS


Seat widths	16", 18" & 20"
Seat depths	16" & 18"
Total weight	15.5 kg.
Weight capacity	120 kg. / 265 lbs.
Casters	8" x 2"
Rear wheels	24" pneumatic, optional Solid Wheels Quick-release
Armrests	Height adjustable
Footrests	Removable & adjustable
Backrest	Tension adjustable
Axle	Height and COG adjustable
HCPCS	K0003

Folding Frame



High Quality Cross

Quick Release  
Rear Wheels



**Shipping dimensions**  
A 33"  
B 13"  
C 39"

**Weight 16 kg.**

## BENEFITS

For people who require a light and safe transfer chair.

## Discover our products **ONE**



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product**  
**Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)

