

Cushion FUSION

JAY®



Optimized for superior skin protection while maximizing stability and comfort. It combines the skin protection capabilities of the J2 with the maximum stability of the J3.

It can be ordered with the JAY Flow fluid, CryoFluid, or air flotation inserts, which adjust to each individual's shape. **Image 1.**

CryoFluid actively cools the settled skin surface while evenly distributing pressure, reducing shear and reducing the risk of wetness. **Image 2.**

Multi Layer, with Contoured Foam base. **Image 3.**

Double deck system. The inner cover is waterproof and easy to clean. **Image 4.**



Image 1



Image 2



Image 3



Image 4




For you. The best technology.

Wheelchair cushions designed to provide the user with support, comfort and protection.

SPECIFICATIONS

Width range	14" to 24" / CRYO 16" to 22"
Length	14" to 24" / CRYO 16" to 22"
Product weight	2.4 kg. / 5.3 lbs. (16x18)
Height	10.6 cm / 4" (3" 3/8" *RP)
Weight capacity	14"-21" 136 kg. / 300 lbs. 22"-24" 226 kg. / 500 lbs.
Inner cover	Moisture resistant AquaGuard®
Outer cover	Stretch™
Settings	In sizes
Options	Fusion / Core / *Reduce Profile

Three options
in height
and accessories



Optimum pelvic
immersion area

Belcros
to keep
the inserts



BENEFITS

For active users with lack of sensitivity and stability

Discover our products **ONE**



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com

