

Compact Scooter

limosa 



 **ONE**

Compact enough to travel with you but rugged enough to take outdoors, Limosa® scooter combines the performance of a sport scooter with the portability of a travel scooter.

Limosa® combines larger wheels, higher ground clearance, more power and a top speed. **Image 1.**

Detachable 5-piece design fits in most trunks or rear seats. **Image 2.**

Height adjustable seat and rudder adjustable ensure you will have a comfortable ride. **Image 3.**

The external charger makes charging easy. Directly connect or carry box of the battery to a place to charge. **Image 4.**



Image 1



Image 2




Image 3



Image 4



Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

Frame	Detachable 5 pieces
Seat Width	17"
Seat Deep	15"
Backrest height	15"
Product weight	46.5 kg. / 102.5 lbs.
Weight capacity	136 kg. / 300 lbs.
Wheels	8"x 2" Solid
Base to floor height	4"
Armrests	Folding and padded
Speed	6 mph
Color	Green oil



Shipping dimensions

A 47"

B 20"

C 21"

Weight 50 kg.

BENEFITS

Designed to be driven on sidewalks and pedestrian areas.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com

