

# QRi

 QUICKIE®



The QRi® is an ultra lightweight rigid wheelchair with a strong open frame design perfect for your on-the-go lifestyle. Designed for form and function, the QRi is a hybrid of durability and low-maintenance style at an unbelievable value.

Fluid caster design, and optional WC-19 approved transit. **Image 1.**

The QRi wheelchair comes standard with 0" or 1" of frame inset. **Image 2.**

Including the Folding Lock-Down Angle Adjustable Back. **Image 3.**

QR Quick extraction system on rear wheels completely safe for better transportation and assembly. **Image 4.**



Image 1



Image 2



Image 3





Image 4

# For you. The best technology.

For those who require cutting-edge technology to keep up with the demands of their active life.

## SPECIFICATIONS

Frame	Aluminum 6062
Seat Width	12" to 20"
Seat Depth	12" to 20"
Product weight	10 kg. / 22 lbs.
User Weight Capacity	120 kg. / 265 lbs.
Gravity center	Adjustable 0.5" up to 3"
Accessories	Optional
Settings	Optional
Available colours:	



Compatible with  
JAY cushions and  
backrests



## BENEFITS

For active people mainly spinal and motor injuries.

## Discover our products



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

