

WAVE

Bath system



When we make bathing easier, our clients are healthier. The new WAVE, a bath chair that improves the bathing and showering experience for both the client and the caregiver.

Durable and soft to the touch material, it has antimicrobial properties, is waterproof and latex free. **Image 1.**

Adjust the angle of the backrest and seat up to 75° with one hand, for greater safety and peace of mind. **Image 2.**

The shower support allows the user to be placed in the shower at a height that allows for better cleaning and protects their back. Locking wheels (4") allow for easy, smooth mobility. **Image 3.**

The ingenious tub transfer base is ideal for large or heavy users who need a shower over the tub. **Image 4.**



Image 1



Image 2



Image 3



Image 4




Adapted equipment. **Synonym of service.**

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS

Total user measurement	76 cm - 188 cm 30" - 74"
Frame measurements	Small, Medium & Large
Seat width	15"/15"/17.5"
Backrest height	23"/26.5"/32"
Max user weight	34 kg. - 113 kg. 75 lbs. - 250 lbs.
Product weight	4.6 kg. - 14.5 kg. 10.2 lbs. - 32 lbs.
Use	Transfer bath system
Accessories	Optional
Settings	No tools required
Fabric colors	

Required user
measurement
(Total height)



Bathtub support
raises the base
up to 30cm.

Easy to store
and transport



BENEFITS

Suitable for patients with neuronal and psychomotor damage or cognitive psychic disability.

Discover our products **ONE**



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

