

# Reclining chair



**coyote<sup>R</sup>**

**ONE**

The OneLoh Coyote R is a reclining, folding wheelchair at a great price. It offers a quick-release axle, removable armrests, and elevating legrests.

Folding frame for easy transport and storage with locking bar for rigidity during recline. **Image 1.**

Legrests elevate from 50° to 0° while also being height adjustable and removable for transfers. **Image 2.**

Reclining backrest from 85° to 175°. **Image 3.**

Quick release rear wheels for easy transportation. **Image 4.**




Image 1




Image 2



Image 3



Image 4




# Keeping you in mind. More options.

Quality products with the required accessories at a great price.

## SPECIFICATIONS

Frame structure	Steel Silver vein
Seat width	16", 18" and 20"
Seat depth	16"
Product weight	26 kg
Weight capacity	136 kg / 300 lbs
Front wheels	8"x2"
Rear wheels	24" pneumatic quick extraction
Armrest	Detachable
Footrest	Elevatable, foldable and removable
Backrest height	20" / With head support 31"
Head support	Including



## BENEFITS

A special needs chair that provides recline and positioning.

## Discover our products



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product**  
**Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

