

Wheelchair with 45° tilt

 Coyote T

 ONE



The Coyote T wheelchair is a manual, aluminum wheelchair that features a manual tilt of 45°. It is designed to provide users with a comfortable and customizable seating position, making it an ideal choice for individuals who require additional support or pressure relief. The chair's lightweight, aluminum frame makes it easy to maneuver and transport. Its manual tilt-in-space system allows for precise adjustments to the seating position, providing maximum comfort and support.

Folding frame for easy transport and storage. **Image 1.**

Removable abductor for improved positioning. **Image 2.**

Red and blue color options and optional tray table. **Image 3.**

Tilt is controlled by levers on the attendant push bars. **Image 4.**



Image 1



Image 2



Image 3



Image 4

Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

Frame	Aluminum, foldable
Seat width	14", 16"
Seat depth	14"
Backrest height	18" / 22" with headrest
Product weight	21 kg
Weight capacity	90 kg / 198 lbs
Front wheel	6" solid
Rear wheel	16" solid
Seat, backrest	Breathable padding
Armrest	Adjustable in height
Foot support	Elevatable and removable
Head support	Padded



Shipping Measurements

A 32"

B 16"

C 27"

Weight 27 kg.

BENEFITS

A special needs chair that provides tilting and positioning support.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

