

Folding Power Chair

FELIX



ONE

The OneLoh Felix folding power chair is the perfect blend of power, portability, and style. The dual motors are controlled by a PG controller and powered by a 2x12 V battery box that can be easily removed for folding, storage, and transit.

Folding frame for easy portability and storage for users with limited space. **Image 1.**

The easy-to-use 4-button PG joystick controls dual Linux engines. **Image 2.**

Removable battery box. **Image 3.**

Disassembles for easy transport. **Image 4.**




Figure 1



Figure 2



Figure 3



Figure 4



Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

Frame	Steel, Folding
Seat Width	17", 18" & 19"
Seat Deep	17"
Backrest height	16"
Backrest type	Tension adjustable
Product weight	65 kg
Weight capacity	136 kg / 300 lbs
Wheels	Rear 12" solid. Casters 8"x 2" solid.
Armrests	Padded
Speed	4 mph / 7 kmh
Battery	Lead Acid 12V/25AH sealed
Color	Blue
HCPCS	K0800 & K0815



Shipping dimension

A 34"

B 15"

C 39"

Weight 62 kg.

BENEFITS For users with reduced mobility partial or total lower limbs.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

