

OneLoh Cushion



The OneLoh Basic Cushion offers an ergonomic design that adapts to the contours of the body, making it more comfortable due to improved positioning of the pelvis. The waterproof, stretch, anti-friction cover's invisible zipper prevents scratches and facilitates changing or washing the cover.

Nylon cover. **Image 1.**

High rebound contoured foam. **Image 2.**

Ventilation holes. **Image 3.**

Soft and breathable cover. **Image 4.**




Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

Measures	14x14", 16x16", 18x16", 18x17", 20x18"
Seat thickness	8 cm / 3.1"
Lower thigh thickness	7 cm / 2.7"
Medial thigh thickness	4.5 cm / 1.7"
Lateral thigh support	4.5 cm / 1.7"
Pelvic Contour	3.5 cm / 1.3"
Weight capacity	130 kg. 286 lbs.



Shipping dimensions

A 17

B 3.2"

C 18.2"

Weight 5 kg.

BENEFITS

Auxiliary for positioning and user comfort.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

