

Q300 M Mini



The Q300 M Mini is a narrower true center-drive motorized wheelchair, with an ultra-small 17" turning radius. All-wheel suspension dampens each wheel independently. Overcome obstacles up to 3" with drive wheels of 12".

It offers the recommended option for JAY cushions and backrests tailored to the user's specific needs. **Image 1.**

LUCI Smart Technology prevents tips, collisions, and falls. Optional add-on for SEDEO PRO seating. **Image 2.**

C-Me® allows you to rise up to 12", while traveling up to 4.8 kph. **Image 3.**

Allows for indoor and outdoor use thanks to the base accessories, front and rear wheels. **Image 4.**




Image 1



Image 2



Image 3



Image 4



For you. The best technology.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Seat width	10" a 16"
Seat depth	10" a 19"
Backrest height	12" a 20"
Weight capacity	75 kg / 165 lbs
Product weight	From 114 kg
Wheel system	Spyder Track 12"/14" and 7"
Features to choose	Tilting, Lifting, C-Me Recline, Elevation of feet
Motors	4 poles
Speed	6 mph standard
Bluetooth	Optional

Available colours:



BENEFITS

For active people mainly spinal and motor injuries.

Electronic functions
that can be configured
and customized



2 sealed batteries

SEDEO PRO
Configurable Seating
with JAY Syncra



33" short base length

Discover our products



DESIGNED FOR LATIN AMERICA

Our products are designed to meet the needs of Latin America and the Caribbean, offering quality and good prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and motorized wheelchairs, to mobility scooters and standing chairs

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com

