

Soft-Touch Liners



They offer additional support and comfort and are recommended for children with mild to moderate sitting support needs. It is made of a soft-touch material that is latex-free and antimicrobial. The seat and back cushions can be used together or separately.

Soft contoured seat cushion and/or back cushion. **Image 1.**

They can be attached to most seating systems. **Image 2.**

It has been designed and molded with a "hinge" to move comfortably and support the body while reclining. **Image 3.**

They are easily transported. **Image 4.**

www.lohmedical.com




Image 1



Image 2




Image 3



Image 4



Support. Solutions for you.

Clinically sound, innovative design and exceptional value.

SPECIFICATIONS

Sizes	Size 1 to 5
Seat cushion width	20 - 40.6 cm / 8 - 16"
Deep seat cushion	20 - 40.6 cm / 8 - 16"
Backrest width	18.75 - 35.5 cm / 7.5 - 14"
High backrest	25 - 40.6 cm / 10 - 16"
User weight	9 kg. up to 91 kg. / 20 - 200 lbs.
Colors	Cherry Dark gray Lilac Aqua

3 Point Pelvic
Cushion Harness
and Trunk Harness.



Built-in abductor
for seat cushion

For special recliner
support system



BENEFITS

They are intended for children with mild to moderate positioning needs.

Conheça nossos produtos



PROJETADOS PARA NOSSOS TERRITÓRIOS

Nossos produtos são projetados para atender às necessidades de cada um de nossos territórios, oferecendo qualidade e preços acessíveis.



VERSÁTIL E FUNCIONAL

Dispomos de uma variedade de produtos que vão desde uma vasta gama de cadeiras de rodas manuais, motorizadas e scooter de mobilidade.

Para mais informações sobre como obter este produto
Entre em contato com seu representante local da Loh Medical

Encontre seu representante local em www.lohmedical.com

