

ANCHOR Bath chair



Provides a worry-free bathing experience for children and caregivers alike. We make safety and comfort a priority. You can make angle, height and depth adjustments allowing for long term growth and use.

Folds flat for easy storage and transport. **Image 1.**

Angle adjustable back, seat and calf support can be configured for individual requirement. **Image 2.**

Positioning straps and head support can be placed anywhere a long length of the chair. **Image 3.**

Optional accesories. **Image 4.**



Image 1



Image 2



Image 3




Image 4



Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS	
Seat width	XS - 12" S - 14.5" L - 16"
deep seat	XS - 10.5 - 13" S - 10.5-13" L - 13-16"
Backrest height	XS - 18.5" S - 24" L - 32"
Backrest width	XS - 14.5" S - 17" L - 18.5"
foot support	XS - 10.2" S - 11" L - 14"
Backrest recline	0°, 30°, 60°, 90°
Angle support feet	30°, 60°, 90°
Max user weight	70 kg / 155 lbs
Product weight	XS - 5.8 kg / 12.75 lbs S - 6.6 kg / 15.5 lbs L - 8.2 kg / 18 lbs
Fabric color	Blue



BENEFITS

Suitable for patients with neuronal and psychomotor damage or cognitive psychic disability.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

