

Adaptive Desk



This desk is the inclusive learning solution that adjusts to every student's needs in seconds. With a wide, wheelchair-friendly frame and a large range of height adjustment, it accommodates all learners with ease.

Height adjustment: 23"–43". **Image 1.**

Optional 0°–50° surface tilt provides customization. **Image 2.**

The casters provide mobility and can be locked for stability. **Image 3.**

Anti-tip rollers keep the desk from tipping forward. **Image 4.**



Image 1



Image 2



Image 3




Image 4



Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS	
ITEM	R250
Height	23" to 43" / 58 to 109 cm
Frame Length	36.5" / 93 cm
Frame Depth	28" / 71 cm
Product Weight	68 lbs / 31 kg
Maximum Load	110 lbs / 50 kg 55 lbs / 25 kg when tilted
Elbow Supports	Supports 25 lbs / 11 kg
Adjustments	No tools required
Table Colors	  



BENEFITS

Desk for students with special needs.

Discover our products



DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

