

Allstar

RGK



The AllStar's simple adjustability teamed with the omission of unnecessary additional weight provides the perfect balance of light weight and flexibility. It's the link between the beginner's Club Sport and the more professional Elite. Player development requires something that can be fine-tuned but also still made to your own specific measurements.

Frame Options. Single or double anti-tip. Fixed or height adjustable backrest. **Image 1.**

Backrest center of gravity adjustment system. **Image 2.**

Upholstery adjustable tension straps for better support and positioning. **Image 3.**

Wide range of ratchet and Velcro strapping to connect you securely to the frame. **Image 4.**



Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Frame	7020 Aluminum
Frame Type	Rigid
Seat Width	10" to 20" / 25 cm to 50 cm
Seat Depth	10" to 20" / 25 cm to 50 cm
Center of Gravity	Adjustable from 0" to 7"
Weight Capacity	100 kg / 220 lbs
Product Weight	From 9.8 kg / 21.6 lbs
Accessories	Optional
Adjustments	In measurements
Available Colors:	27 options



BENEFITS

For active people mainly motor injuries.

Discover our products



DESIGNED FOR OUR TERRITORIES



Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

