

Moldable seat Beanseat



This large orthopedic seat stabilizes the pelvis, torso and head. It enables the patient to be placed in a deep sitting or semi-lying position, stabilizing the whole body.

The materials are certified non-flammable, antibacterial and latex-free.
Image 1.

Two sizes of width, height and length.
Image 2.

Includes manual air pump, neoprene tape and other optional accessories.
Image 3.

Flexible cover that increases the level of hygiene. It is impermeable to liquids and permeable to water vapor and air.
Image 4.




Image 1




Image 2



Image 3



Image 4



Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS	
Width	Medium 60 cm Large 70 cm
Length	Medium 80 cm Large 100 cm
Rear height	Medium 50 cm Large 70 cm
Front height	Medium 40 cm Large 50 cm
User Weight capacity	130 kg / 286.6 lbs
Product weight	10 kg / 22 lbs
Accessories	Optional



Pelvic stabilizer



Transfer platform

BENEFITS

Designed for users who require moderate positioning and stability.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

