

Buzzaround LT Scooter

 GOLDEN



It is intended for use by people of all ages who may have difficulty walking distances or for periods of time.

This sporty design features a lightweight stadium-style seat and a 300 lbs. weight capacity. **Image 1.**

All Golden scooters can be used with left or right controls and swivel operation. **Image 2.**

FDA approved class II medical devices. **Image 3.**

Optional accessories for the user and for the scooter. **Image 4.**



Image 1



Image 2



Image 3



Image 4

Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

ITEM	GB107
Seat width	17"
deep seat	16"
Backrest height	16"
Product weight	45 kg / 98 lbs
Weight capacity	135 kg / 300 lbs
Wheels	8" x 1.5" Solid
Base to floor height	6 cm / 2.5"
Speed	4 mph
Color	Red



BENEFITS

Designed to be driven on sidewalks and pedestrian areas.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

