

CHROME

RGK



CHROME is a perfect example of how reducing excess unnecessary parts and understanding the materials you are working with can allow for the creation of a great chair. This chair is an incredibly lightweight and durable model, and will do wonders for your daily needs.

Chrome moly steel frame, 25mm diameter frame tube size. **Image 1.**

Choice of open or closed frame design. **Image 2.**

Adjustable upholstery with adjustable tension straps for better support and positioning. **Image 3.**

Fixed or adjustable footrest position. **Image 4.**



Image 1



Image 2



Image 3



Image 4


Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

| | |
|----------------------|---------------------------------------|
| Frame Material | Chrome moly steel frame |
| Frame Style | Rigid |
| User Weight Capacity | 125 kg / 275.5 lbs |
| Gravity center | Adjustable or Fixed |
| Product Weight | From 9.1 kg / 20 lbs |
| Seat Width | Bariatric up to 21.5" (55 cm) |
| Caster | 4" or 5" with aluminum fork |
| Backrest | Fixed or folding adjustable in height |

Add personal embroidery on the backing.



Rear wheels
Quick release

BENEFITS

For active people mainly motor injuries.

Discover our products **ONE**



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

