

# COMPASS

 rifton.



It is the perfect inclusion chair. With its supportive seat, sturdy armrests, and stable frame, it offers the support they need to feel safe and comfortable.

Standard legs have two inches of height adjustment. **Image 1.**

Pelvic harness secures the hips and upper thighs without pressure on the abdomen. Provides a secure and comfortable base for correct posture. **Image 2.**

High armrests provide support and aid in transfers. **Image 3.**

Stability feet widen the base of support and prevent tipping. 2" height adjustment. **Image 4.**



Image 1



Image 2



Image 3






Image 4



# Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS	
User measurement	91 cm to 160 cm / 36" to 63"
Frame Sizes	1 - 5
Seat width	27 cm to 46 cm / 10.5" to 18.2"
Seat depth	24 cm to 42 cm / 9.5" to 16.5"
Max user weight	113 kg / 250 lbs
Use	Inside
Accessories	Optional
Settings	No tools required
Color	Blue

Required user measurement (total height and knee to heel)	
	Pelvic belt for safety and posture
	
5 sizes available	

## BENEFITS

Suitable for patients with neuronal and psychomotor damage or cognitive psychic disability.

## Discover our products



### DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product**  
**Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

