

Club Sport

RGK



The sports wheelchair that allows you to get active without a major commitment. Only the best materials are used in your Club Sport. Chrome Moly is famous for its strength, durability and is surprisingly not that heavy when worked with properly.

Configurable sizes. **Image 1.**

Suitable for multisport use - e.g. Basketball, Tennis, Rugby. **Image 2.**

Adjustable height backrest and centre of gravity. **Image 3.**

Quick release single anti-tip. **Image 4.**



Image 1



Image 2



Image 3




Image 4



Keeping you in mind. More options.

Quality products with the required accessories at a great price.


SPECIFICATIONS

Seat width	13" to 17"
Seat depth	12" to 17"
Frame height	19" to 21" (front)
Frame height	16" to 21" (rear)
Product weight	From 10 kg / 22 lbs
Weight capacity	125 kg / 275 lbs
Gravity center	Adjustable
Camber	18°
Casters	3"
Accessories	Optional
Available colours	23 available

Regulatory measures



Rear wheels
Quick release



BENEFITS

For active people mainly motor injuries.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

