

Comfortable Backrest



The Comfortable stabilizing backrest provides comfortable support and relief at the same time. The modifiable shape offers personalized adjustment to the user's needs, correcting postural defects and thus improving the control of the user's body.

The materials are certified non-flammable, antibacterial and latex-free.
Image 1.

Three sizes of width and height.
Image 2.

Includes manual air pump, and other optional accessories.
Image 3.

The anatomical shape of the backrest and light weight allow for easy installation in a wheelchair, rehabilitation seat, chair or car.
Image 4.




Image 1




Image 2



Image 3



Image 4



Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Width	Small 52 cm
	Medium 56 cm
	Large 70 cm
Height	Small 40 cm
	Medium 50 cm
	Large 50 cm
User Weight capacity	100 kg / 200 lbs
Product weight	1.5 kg / 3.3 lbs
Accessories	Optional

Flexible tube for greater stability and rigidity



BENEFITS

Designed for users who require moderate positioning and stability.

Discover our products



DESIGNED FOR OUR TERRITORIES



Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

