

Aluminum Wheelchair



 **ONE**

Designed for comfort and mobility, this lightweight folding aluminum wheelchair features adjustable components and quick-release wheels for easy transport and compact storage. Ideal for: Patients with any mobility need who value flexibility and daily independence.

Folding aluminum frame for easy transfer and storage if space is limited.

Image 1.

The armrest are height adjustable to accommodate all users. **Image 2.**

The axle can be adjustable in two dimensions: height and depth (center of gravity). **Image 3.**

The footrest hangers are adjustable in height and can be easily removed. **Image 4.**

Image 4.



Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS


Seat widths	16", 18" & 20"
Seat depths	16" & 18"
Total weight	15.5 kg
Weight capacity	120 kg / 265 lbs
Casters	8" x 2"
Rear wheels	24" pneumatic, optional Solid Wheels Quick-release
Armrests	Height adjustable
Footrests	Removable & adjustable
Backrest	Tension adjustable
Axle	Height and COG adjustable
HCPCS	K0003

Folding Frame



High Quality Cross

Quick Release
Rear Wheels



Shipping dimensions
A 33"
B 13"
C 39"

Weight 18 kg.

BENEFITS

For people who require a light and safe transfer chair.

Discover our products **ONE**



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com



ver_2025.11.11E

