

# Tilt Wheelchair 45°



This pediatric aluminum wheelchair includes tilt and recline functions, with postural components and adjustable features that support growth. Ideal for Children needing daily mobility support with a safe design.

Folding frame for easy transport and storage. **Image 1.**

It features a removable abductor for improved positioning. **Image 2.**

The frame is available in red and blue, and tray. **Image 3.**

The tilt is activated by a completely safe lever on the handlebar. **Image 4.**



Image 1



Image 2



Image 3




Image 4

# Keeping you in mind. More options.

Quality products with the required accessories at a great price.

## SPECIFICATIONS

Frame	Aluminum, foldable
Seat width	14", 16"
Seat depth	14"
Backrest height	18" / 22" with headrest
Product weight	21 kg
Weight capacity	90 kg / 198 lbs
Front wheel	6" solid
Rear wheel	16" solid
Seat, backrest	Breathable padding
Armrest	Adjustable in height
Foot support	Elevatable and removable
Head support	Padded



### Shipping Measurements

A 32"

B 16"

C 27"

Weight 25 kg.

## BENEFITS

A special needs chair that provides tilting and positioning support.

## Discover our products



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

