

Wheelchair with Manual Tilt and Recline



We now offer a rigid frame, reducible chair that has tilting up to 45° and also manual recline included, high quality, at affordable prices, custom configured as required for children and adults.

Tilt from -3° to 45°. **Image 1.**

Recline from 90° to 125°. **Image 2.**

Solid 16" rear wheels, with drum brake. **Image 3.**

Leg supports adjustable in angle, height and tibiotarsus. **Image 4.**

Other features:

- Single angle adjustable handlebar.
- Easy-to-use retractable anti-tip
- Desk-type armrests, adjustable in height and removable.
- Aluminum seat panel and adjustable hip supports.
- Includes transit system.
- The backrest and cushion are not included in the basic configuration, any positioning system of your choice can be attached, consult our team who will advise you on the best cushion and backrest to attach.




Image 1




Image 2



Image 3




Image 4

Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

Frame	Aluminum, Rigid reducible Black color
Seat width	14" up to 19"
Seat depth	14" up to 19" (Aluminum panel)
Product weight	34 kg
Weight capacity	136 kg
Front wheel	7" solid
Rear wheel	16" solid
Armrest	Adjustable in height and removable
Head support	Anatomical, adjustable in angle and depth
Backrest and cushion	Not included
Trunk supports	Not included



Shipping Measurements

A 30"

B 32"

C 33"

Weight 40 kg



BENEFITS

A wheelchair designed to suit special needs that offers tilting, reclining and that adapts to all ages

Discover our products ONE



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com



ver_2025.11.11E

