

# Cushion GS



The JAY® GS® Wheelchair Cushion incorporates a highly durable, stable foam base and JAY's proven fluid insertion technology with the unique ability to grow with the client while maintaining adequate support, making it an excellent solution for growing children and young adults.

Contoured foam base with built-in growth and easy modifications.  
**Image 1.**

Available with a wide selection of 2" and 3" positioning components.  
**Image 2.**

It has a FLUID insert that requires no maintenance and adapts to the shape of each individual.  
**Image 3.**

Stretch cover with Ionic+™ fiber and 3DX spacer fabric.  
**Image 4.**



Image 1



Image 2



Image 3



Image 4




# Keeping you in mind. More options.

Wheelchair cushions designed to provide the user with support, comfort and protection.

## SPECIFICATIONS

Width range	10" a 18"
Depth	10" a 20"
Product weight	3.1 kg / 6.9 lbs
Height	10 cm / 4" with components
Weight capacity	14"-21" 136 kg / 300 lbs
Foam base	Precontoured foam
Insert	Flow™ fluid compartmentalized pad
Outer cover	X-Static Stretch™
Settings	Optional
Options	2" and 3" components

Belcros to maintain and position components



Optimum pelvic immersion area

## BENEFITS

It is an excellent positioning solution for growing children and young adults.

## Discover our products **ONE**



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

