

Cushion X2

JAY®



The JAY X2 is a lightweight, maintenance-free, fluid wheelchair cushion with a reduced profile that is designed to provide stability and pressure relief to active clients. Using elements from the JAY Xtreme and Active cushions, the X2 combines JAY's time-tested designs with its latest technologies.

Light and stable foam base with reduced profile. **Image 1.**

It features a small JAY Flow fluid insert, which adapts to each individual's shape and secures bony prominences to help protect the skin from injury. **Image 2.**

Innovative double deck system. **Image 3.**

The bag in the photo is optional. **Image 4.**



Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

Wheelchair cushions designed to provide the user with support, comfort and protection.

SPECIFICATIONS

Width range	14" a 20"
Depth	14" a 20"
Product weight	1.7 kg / 3.8 lbs (16"x16")
Height	2.72 cm / 2.5"
Weight capacity	136 kg / 300 lbs
Foam base	Foam
Insert	JAY Flow™ Fluidio Tripad
Inner cover	Waterproof with anti-wick seam
Outer cover	Stretch™
Settings	In measurements

Replacement fluid for each cushion size.



The inner cover is waterproof reinforced with anti-absorption thread

BENEFITS

For active users with lack of sensitivity and stability

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

