

# EZ-2 Pool Lift



The total solution for access to pools and spas. For residential, commercial and medical use.

Anchor: Choose from concrete or wood deck applications. **Image 1.**

Cover: Protects your investment from the elements. **Image 2.**

4 Point Suspension Option – Safer and able to be used with many standard lanyards. **Image 3.**

Hard seat – additional seating option for a wider range of users. **Image 4.**



Image 1



Image 2



Image 3



Image 4



# Keeping you in mind. More options.

The leaders in recreation, fitness and accessibility, safe and guaranteed.

## SPECIFICATIONS

ITEM	F-EZPL2
Seat range	31 - 50 cm / 12" - 20"
Weight capacity	181 kg / 400 lbs
Product weight	31.7 kg / 69.9 lbs
Rotation range	360°
System	Hydraulic manual
Accessories	Stainless steel chain Seat models (2) Anchor (not included) F-033EZA – Standard anchor



## BENEFITS

Manufactured for pools and spas, they are ideal for both commercial and residential use.

## Discover our products



### DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

