

Easy Cushion



The JAY® Easy™ is a great value wheelchair cushion with JAY Flow fluid for skin protection and positioning. The Easy™ cushion features a hi-resiliency, contoured foam base that accommodates a curved or flat seating surface and JAY Flow™ fluid tri-pad.

Contoured Foam Base with Optional Curved Bottom. **Image 1.**

It has a small fluid insert, which adjusts to the shape of each individual. **Image 2.**

VISCO base option. **Image 3.**

Can be ordered with colored outer shell side bands. **Image 4.**



Image 1



Image 2



Image 3



Image 4

Keeping you in mind. More options.

Wheelchair cushions designed to provide the user with support, comfort and protection.

SPECIFICATIONS

Width range	14" to 24"
Length	14" to 24"
Product weight	1.8 kg / 4 lbs (16"x18")
Height	7.6 cm / 3"
Weight capacity	113 kg / 250 lbs
Foam base	Pre-shaped foam
Insert	VISCO fluid or model
Outer cover	Microclimate / Incontinence
Settings	In width and depth

BENEFITS

Patient with symmetrical or mild asymmetrical posture requiring moderate postural support

Discover our products



DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

