

Extreme X8



magic
Mobility

4x4 Standar Powered Wheelchair. Unrivalled off-road performance and reliability make it perfect for outdoor adventurers.

With a 12" seat lift, the C-Me electric lift function. (5) Travel at 5 km/h with 12" seat elevation. **Image 1.**

Four independent, high speed, high torque and high efficiency motors. **Image 2.**

Hinged front frame with low-pressure knobby tires. **Image 3.**

Offers five standard colours and three specialised prints. **Image 4.**



Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Seat width	13" up to 24"
Seat depth	13" up to 24"
Product weight	From 145 kg
Weight capacity	154 kg
Backrest height	16.5" up to 24"
Maximum rank	33 km
Maximum speed	10 kmh
Wheels	14"
Turning radius	66 cm
Total width	28"



- Seat lift and tilt 12" / 45°
- 95-180° Legrest elevation
- Tilt 50°
- 95-160° Recline

BENEFITS

For users with partial or total reduced mobility of lower limbs.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

