


Floating beach Wheelchair



The Floating Beach Wheelchair is a lightweight, sturdy, all-terrain wheelchair that provides a seamless transition from boardwalk-to-beach-to-water. The wheelchair is designed for access in sand, soft soils and snow. Durable and easy to maintain, this floating wheelchair is constructed of lightweight stainless steel with an aluminum frame, mildew resistant mesh seat and puncture resistant balloon tires.



Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS

Wheels	12" and 16" pneumatic tires
Seat width	18"
Seat depth	17"
Backrest height	17"-19"
Weight capacity	115 kg / 250 lbs
Product weight	25 kg / 55 lbs
Armrests	Floats
Wheel material	Yellow (RUBBER)
Hardware	Stainless steel
Accessory	Safety harness



BENEFITS

Excellent help for patients with different transfer needs and limited mobility.

Discover our products **ONE**



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

